

Bacillol® Wipes

Alcohol-based disinfection wipes



Aldehyde-free, alcohol-based disinfection wipes for rapid disinfection of alcohol-resistant surfaces. Developed for the use in combination with BODE X-Wipes dispenser system.

Bacillo® Wipes

Product properties

- pre-soaked disinfection wipes for the BODE X-Wipes dispenser system
- broad spectrum of effect
- short action time within 30 seconds (VAH)
- compatible with a wide range of alcohol-resistant surface materials
- good coverage and rapid drying
- tear-proof, absorbent and lint-free fleece
- use-solution usable for 4 weeks after opening

With Bacillo® Wipes used in combination with the BODE X-Wipes dispenser system the user will benefit from many advantages: the one-time use of the wipes reduces the preparation time and guarantees a high hygiene standard. At the same time, the high-quality materials provide for low consumption of fleece wipes and disinfectant. With the dispensing system properly closed, the working solution can be used for 4 weeks after opening.

Composition

Active ingredients of the solution: Propan-1-ol 450 mg/g, propan-2-ol 250 mg/g, ethanol 47 mg/g.

Microbiological activities

Bactericidal, yeasticidal, fungicidal, tuberculocidal, mycobactericidal, virucidal against enveloped viruses (incl. HBV, HIV, HCV), adenoviruses, FCV, MNV, polyoma- and rotaviruses.

Areas of application

Bacillo® Wipes used in combination with the BODE X-Wipes dispenser system are suitable for the uncomplicated cleaning and disinfection of:

medical equipment as defined by the Medical Device Directive and other alcohol-resistant parts of medical devices

surfaces that are not sensitive to alcohol according to BPD, such as:

surfaces in hospitals, medical and dental practices, ambulances
surfaces in medical laboratories, residential and nursing homes, etc. and in industrial kitchens and food-processing areas

The use of Bacillo® Wipes in the BODE X-Wipes dispenser system is also advantageous in areas, which require repeated cleaning every day, as the one-time use assures a high level of hygiene as well as a reduced preparation time. In addition, surfaces difficult to access are optimal application areas for Bacillo® Wipes. The wipes can be folded according to the requirements on location and allow a good coverage of even small areas or niches.

Directions for use

Carefully wipe surfaces with Bacillo® Wipes. Ensure complete coverage so that optimum disinfection is guaranteed. The standing time of 28 days (4 weeks) can only be recommended with properly closed protective lid. After withdrawal, the dispensing system has to be closed firmly. Not suitable for the disinfection of invasive medical devices.

Filling of X-Wipes dispenser system with Bacillo® Wipes

1. Open dispenser's lid and insert Bacillo® Wipes sachet.



2. Cut Bacillo® Wipes sachet open with a cutter on the top.



3. Slightly pull out the first wipe from the middle of the roll and draw first Bacillo® Wipe from the roll's middle through the opening of the dispensing system and discard. Put lid on dispenser and press until it snaps in. Afterwards close the dispensing system's protective cover carefully.



4. Complete label and stick it on the indicated area on the X-Wipes dispenser.



Use disinfectants safely. Always read the label and product information before use.



Bacillo[®] Wipes

Dosage

Bacillo[®] Wipes are a ready-to-use product.

Surface disinfection (Wipe on)			
Bactericidal/Yeasticidal ¹⁾			
– dirty conditions	conc.	–	30 sec ²⁾
	conc.	–	5 min
Bactericidal (EN 1040)	conc.	–	1 min
Bactericidal (EN 13727)			
– dirty conditions	conc.	–	15 sec
Fungicidal ¹⁾			
– dirty conditions	conc.	–	5 min
Fungicidal (EN 1275)	conc.	–	5 min
Tuberculocidal ¹⁾			
– dirty conditions	conc.	–	30 sec ²⁾
	conc.	–	5 min
Mycobactericidal ¹⁾			
– dirty conditions	conc.	–	30 sec ²⁾
	conc.	–	5 min
Virucidal against enveloped viruses (incl. HBV, HIV, HCV)			
	conc.	–	30 sec
Adenovirus	conc.	–	1 min
FCV ³⁾	conc.	–	5 min
MNV ⁴⁾			
– clean conditions	conc.	–	1 min
– dirty conditions	conc.	–	1 min
Polyomavirus	conc.	–	10 min
Rotavirus	conc.	–	1 min
Food processing area (Bactericidal+Fungicidal) ⁵⁾ ; Area A			
– less soiled – 20 °C	conc.	–	30 min
– less soiled – 10 °C	conc.	–	30 min
Food processing area (Bactericidal+Fungicidal) ⁵⁾ ; Area A			
– heavily soiled – 20 °C	conc.	–	30 min
– heavily soiled – 10 °C	conc.	–	30 min
Food processing area (Bactericidal+Fungicidal) ⁵⁾ ; Area B			
– less soiled – 20 °C	conc.	–	30 min
– less soiled – 10 °C	conc.	–	30 min
Food processing area (Bactericidal+Fungicidal) ⁵⁾ ; Area B			
– heavily soiled – 20 °C	conc.	–	30 min
– heavily soiled – 10 °C	conc.	–	30 min
Food processing area (Bactericidal acc. to EN 1276);			
– less soiled – 20 °C	conc.	–	1 min
– heavily soiled – 20 °C	conc.	–	1 min
– less soiled – 10 °C	conc.	–	1 min
– heavily soiled – 10 °C	conc.	–	1 min
Fungicidal (EN 1650)	conc.	–	5 min

Listing

All data referring to the active ingredients solution Bacillo[®] AF.

Certification/list issued by the Association for Applied Hygiene (VAH).

CE-labelling in accordance with the Medical Device Directive (MDD).

List issued by the German Veterinary Medical Society (DVG) for Food Hygiene.

Chemical-physical data

All data referring to the active ingredients solution Bacillo[®] AF:

Flash point (DIN 51755)	25 °C
Density (20 °C)	approx. 0.86 g/cm ³
Refraction n _D ²⁰	approx. 1.38

All data referring to the wipe:

Dimensions: 380 mm x 200 mm

Presentation

Sachet with 90 wipes

Labelling in accordance with Directive 1999/45/EC: Xi Irritant. Flammable.

Irritating to eyes. When used correctly, the basic measures of occupational health and safety in accordance with protection level 2/S9 Hazardous Substances Ordinance (GefStoffV) should be observed. At temperatures higher than its flash point of 25 °C Bacillo[®] AF can develop potentially explosive atmospheres.

Keep out of the reach of children. Avoid contact with eyes. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.

The recommendations regarding our preparations are based on scientific tests and are given in good faith. More detailed recommendations, e.g. regarding material compatibility, are only possible in particular cases. Our recommendations are without obligation and do not constitute a warranty. They do not preclude a company's own testing for the intended purposes and processes. In this respect we cannot accept any liability. This complies with our general conditions of sale and supply.

Effectiveness tested with solution

- 1) according to new standard methods; with mechanical action
- 2) rapid disinfection
- 3) surrogate virus for human norovirus, without organic load
- 4) surrogate virus for human norovirus
- 5) DVG

