Gebrauchsanweisung
User manual
Istruzioni per l'uso
Instructions d'utilisation
Instrucciones de uso
Manual de instruc;ões
Navod R pouzítí

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User manual

Please read these instructions carefully before using the device for the first time!

In the event of cardiac arrest, survival and the prevention of long-term neurological damage significantly depend on immediate initiation of resuscitation measures by first aiders. The Cardio First Angel was developed and designed for people without specialized knowledge to overcome their fear of calling out urgently required resuscitation measures and make an easy-to-use aid for cardiac massage available to them. The device is an entirely mechanical device and does not contain any electronic components, which means that the Cardio First Angel is always ready for immediate use. Failure to follow the instructions and safety precautions in this Manual may result in serious injury or death. IMPORTANT NOTE: For every minute of delay, the chance of survival drops by about 10%. However, it is important to understand that the Cardio First Angel cannot always assure survival. In some cases, the cause of the cardiac arrest cannot be survived no matter how quickly the victim is treated.

I. Introduction

The Cardio First Angel is a handy and uncomplicated device which is used in emergency situations to support laypeople in starting the required resuscitation measures on another person. In emergency situations, it is often a serious problem that people without specialized medical knowledge hesitate to take urgently required first aid measures. There can be various reasons for this: the lack of sufficient knowledge, experience, the spontaneous feeling of helplessness, the fear of doing something wrong or the disinclination to come into contact with the patient's body. However, even if laypeople do initiate the resuscitation measures without hesitation, it is great uncertainty about how strongly and firmly or with what frequency you have to press on the thorax. It is not clear where exactly you have to place your hands for the cardiac massage. The Cardio First Angel frees the user from all these considerations. It is simple to use. Another advantage is that direct body contact is avoided.

II. Description of the Cardio First Angel

The Cardio First Angel has very compact dimensions of 110 mm x 170 mm, only weighs 130 grammes and is stored in a sealed in a transparent red plastic box.

The Cardio First Angel is an entirely mechanical device, functioning without power supply and electronic components. The Cardio First Angel is ready for use at all times and without any further

requirements. Please note that the Cardio First Angel is a disposable product that must not be reused again after use once in an emergency situation!

The central component of the Cardio First Angel is the white stable base plate made of PE plastic. Body-embedding polyurethane foam has been applied onto the bottom side of this base plate in order to ensure the uniform pressure distribution and the absorption of fluids. The special outer flange of the Cardio First Angel facilitates correct positioning on the patient's chest.

On the base plate, a complex arrangement of helical springs has been installed, which transmits the compressive force to the thorax. Applying a force of 410 +/- 30 N, corresponding to a CO₂ pressure < 50 mmHg, results in a clearly audible clicking sound.

The Cardio First Angel is completed by a round, signal red top side, serving as pressure application area when used. On this pressure application area, universally understandable pictograms show how the Cardio First Angel is used so that a brief glance ensures safety of use in an emergency situation. However, it is assumed that these detailed instructions for use have already been read and understood in advance.

III. Mode of Use

When performing cardiac massage, the external pressure must be approx. 410 N, which might seem very high for a person. However, this is the only way to achieve the required pressure depth of 50 mm to 60 mm on the thorax. The Cardio First Angel signals when the correct compressive force of ca. 410 N has been reached by a clearly audible metallic clicking sound. When the pressure is released, an equally loud clicking sound can be heard.

This continuous change between the "applying pressure click" and the "releasing pressure click" serves the user as a guide for the rhythm for the cardiac massage. Using this as a rhythmic intermittent change between actively applying pressure and releasing pressure again, using this as a rhythmic indicator achieves the required pressure frequency of approx. 100 to 120 compressions per minute.

The 5'pecial shape of the Cardio First Angel makes it easy to correctly position the device on the lower third of the patient's breastbone. Thus, the pressure can be transmitted optimally to the cardiac muscle inside the bony thorax.

IV. Application

General procedure to be followed if you have to initiate resuscitation measures on an unresponsive:

- (1) Try to speak to the patient and shake them on their shoulder.
- (2) Call for help loudly, draw the attention of other people to the situation and have an emergency call made by dialling the emergency telephone number 112 (Europe) or 911 (USA).
- (3) Carefully lift the chin and stretch the head backwards in order to open the patient's airways.
- (4) Use your senses of vision, hearing or feeling in order to check whether the patient is breathing.
- (5) If you cannot detect normal breathing, no lifting and lowering of the chest with the cardiac ...**

For this purpose, take the Cardio First Angel in your hand, pull the top cover, which is made of transparent red plastic storage box and take out the Cardio First Angel. Remove the device from the patient's thorax until the chest and back are exposed. Knee down next to the patient and turn yourself towards him at the height of their chest. With the foam side facing downwards, place the Cardio First Angel directly on the skin above the breastbone in such a way that the more pointed end of the Cardio First Angel faces towards the patient's feet and approximately halfway ends at the tip of the breastbone. The full lifting force is then applied at a slightly higher point under the red pressure application area and must not be applied to the lower tip of the breastbone, but correspondingly higher (refer to the drawing).

Place your hands on top of each other on the signal red pressure application area and, with your arms extended, start the compression. For this purpose, you must apply such high pressure to the Cardio First Angel with a jerk that you hear the loud clicking sound. Afterwards, quickly release the pressure and exert full pressure again immediately after you hear the next clicking sound. The recommended consisting of applying pressure and releasing the pressure again must follow one another as fast as you hear the click of the two ticking sounds as a rhythm indicator. When doing so, think of the ticking of a clock. You must reach a frequency of 100 to 130 pressure impulses per minute. Your arms should always remain extended in order to ensure that the compressive force can be generated by your shoulder and the upper part of your body. You must continue this resuscitation until you notice the first clear vital signs or until the paramedic or emergency doctor can carry out further resuscitation measures.

To maintain it up, here are the application steps for performing a cardiac massage using the Cardio First Angel again in brief:

- (1) Position the Cardio First Angel on the lower third of the patient's breastbone.
- (2) Make sure to lay your hands on the pressure application area of the Cardio First Angel in such a way that your arms are extended.
- (3) Enough pressure must be exerted so firmly that you hear the clicking sound.
- (4) The alternating rhythm of applying pressure and releasing pressure is provided by the click of the two clicking sounds.
- (5) Perform the cardiac massage with a pressure frequency of at least 100 compressions per minute until you notice the first clear vital signs or until the emergency doctor can take over from you.

V. Limits

The Cardio First Angel in no way replaces the required professional medical treatment and care provided by paramedics and emergency doctors. Therefore, it is absolutely necessary that help is called for even after the resuscitation measures have been started. The Cardio First Angel is only an aid in order to increase the crucial time until the rescue personnel arrives to take the first resuscitation measures. Under no circumstances can it be guaranteed that the resuscitation measures will succeed, since this ultimately depends on the patient's state of health.

VI. Warnings

Use the Cardio First Angel only according to its instructions. Do not use the Cardio First Angel with care and should not be used on children or on the elderly. Do not use the Cardio First Angel in a clean environment and out of the reach of children. Do not expose it to frost or extreme heat. Under no circumstances should you open or disassemble the Cardio First Angel.

Please check at regular intervals (our recommendation: every 6 months) that the Cardio First Angel is functioning properly. In this context, you should also read these instructions for use again in order to safely master using the Cardio First Angel even in stressful situations.

In order to perform a routine functional test, place the Cardio First Angel on a clean and padded ground (e.g. carpeted floor) and firmly press the Cardio First Angel -as described in the instructions above- several times. You can hear the characteristic clicking sounds when pressing and releasing the device.

If the Cardio First Angel falls down or is otherwise damaged, you must replace it with a new device. Hence, in this case, it can no longer be guaranteed that the springs inside are still sound and in the correct position.

Pay attention to the use-by date which, is specified in the month-year format 00 the white base plate of the Cardio First Angel. After this date has been reached, the Cardio First Angel must be replaced with a new device.

Once the Cardio First Angel has been used for resuscitation measures in an emergency situation, it must under no circumstances be used again. It is a disposable device (single use only)!

You can dispose of the Cardio First Angel together with plastic materials, as it is an entirely mechanical device and does not contain any electronic components.

The Cardio First Angel is not to be used on people younger than 16 years! Special care must be taken when the Cardio First Angel is used on weak or elderly people).

Do not use the Cardio First Angel in the proximity of magnetic resonance imaging (MRI)! When using the device together with a defibrillator, observe the safety instructions provided by the manufacturer!

Do not use the Cardio First Angel on people with regular vital signs!

VII. Risks

A special spring inside the Cardio First Angel is designed in such a way that it clicks. Indicate the right

pressure during the cardiac massage. In individual cases, however, it cannot be excluded that further

complications or injuries/damage are possible (e.g. individual ribs broken) even if the Cardio First Angel is used properly. Although the ease of use was one of the main objectives when developing and designing the Cardio First Angel, it is possible, especially in case of partial or severe shock, that the user feels overwhelmed when handling the Cardio First Angel and makes mistakes.

1. We warrant to the original purchaser that the Cardio First Angel will be free of any defect in material and workmanship according to the terms and conditions of the Limited Warranty (Limited warranty). For purposes of the Limited warranty, the original purchaser is deemed to be the original end user of the product. A purchaser of this Limited Warranty is NONTRANSFERABLE and UNASSIGNABLE.

2. LIMITED WARRANTY OBLIGATION: Exclusive Remedy: THE FOREGOING LIMITED WARRANTY IS IN LIEU OF AND SPECIFICALLY EXCLUDES AND REPLACES ALL OTHER EXPRESSED OR IMPLIED WARRANTIES INCLUDING BUT NOT LIMITED TO THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. Some states do not allow limitations on how long an implied warranty lasts, so this limitation may not apply to you. NO PERSON (INCLUDING ANY AGENT, DEALER, OR REPRESENTATIVE OF US) IS AUTHORIZED TO MAKE ANY REPRESENTATION OR WARRANTY CONCERNING OUR PRODUCTS, EXCEPT TO REFER PURCHASERS TO THIS LIMITED WARRANTY. YOUR EXCLUSIVE REMEDY WITH RESPECT TO ANY AND ALL LOSSES OR DAMAGES RESULTING FROM ANY CAUSE WHATSOEVER SHALL BE AS SPECIFIED ABOVE. WE SHALL IN NO EVENT BE LIABLE FOR ANY SPECIAL, PUNITIVE, INDIRECT, CONSEQUENTIAL OR INCIDENTAL DAMAGES OF ANY KIND, INCLUDING BUT NOT LIMITED TO EXEMPLARY DAMAGES, COMMERCIAL LOSS FROM ANY CAUSE, BUSINESS

INTERRUPTION OF ANY NATURE

VIEW OF THE POSSIBILITIES OF SUCH DAMAGES, HOWEVER OCCASIONED, WHETHER BY NEGLIGENCE OR

OTHERWISE. Some states do not allow the exclusion or limitation of incidental or consequential damages, so the above limitation or exclusion may not apply to you.

3. What this warranty does not cover: This Limited Warranty does not cover defects or damages of any sort resulting from, but not limited to, accident, damage while in transit to our service location, product tampering, unauthorized use, unauthorized service, unauthorized product use, opening, failure to follow instructions, improper use, abuse, neglect, fire, flood, war or acts of God. We make no warranty claim as to the condition of the Cardio First Angel with any other product, part or accessories.

For further information:

The manufacturer reserves the right to make technical changes to the product as a result of the further development. The Cardio First Angel is a class I medical device. OIMD DE/CA64/128442/01, GMDN 60&41

Symbols used on the white base plate of the Cardio First Angel:

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


O

O

• It is important to understand how to use your Cardio First Angel safely. Please read these warnings and carefully.

• FLUIDS. Do not attempt to get into the Cardio First Angel. Avoid spilling any fluids on the Cardio First Angel. Spilling fluids into the Cardio First Angel may damage it. A damaged Cardio First Angel may cause the device to perform improperly, and/or injure the victim or the user.

• CHILDREN. Keep the Cardio First Angel out of reach of children to avoid the potential risk of damage to the device.

• DEVICE HANDLING. The Cardio First Angel was designed to be sturdy and reliable for many different use conditions. However, handling the Cardio First Angel too roughly can damage it.

• MAINTENANCE. Improper maintenance may damage the Cardio First Angel or cause it to function improperly. Maintain the Cardio First Angel according to the specifications as set forth in this Manual.

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Inotech Kunststofftechnik GmbH, Schtschstr. 3, 92504 Nabburg, Germany