

| SKIN CONCERN | | KEY INDICATION | PROGRAMME | PRE-SET PROTOCOLS | | | TIME | JOULES |
|------------------------------------|---|-------------------|-----------|-------------------|--------------|--------------|------------|--------|
| | | | | BLUE 415NM | RED 633NM | NIR 830NM | | |
| SKIN REJUVENATION | Fine lines, photo-damage, loss of elasticity, enlarged pores, uneven skin tone | REJUVENATION | | | H | H | 20 MINUTES | 192 |
| COMPLEXION | Dull and sallow skin, dehydration, lack of radiance, puffiness | REJUVENATION | | | H | H | 20 MINUTES | 192 |
| DRY SKIN | Dry and tight, flaky or rough skin texture, easily irritated, sensitivity and facial redness | INFLAMMATORY SKIN | | | | M | 20 MINUTES | 37.2 |
| PIGMENTATION: PHOTO-DAMAGE | Photo-damage, sun-induced age spots, uneven skin tone | REJUVENATION | | | H | H | 20 MINUTES | 192 |
| PIGMENTATION | Hormonal, post inflammatory hyperpigmentation (PIH), medication, photo-sensitivity | PIGMENTATION | | | H | | 20 MINUTES | 126 |
| ACNE: MILD TO MODERATE | Acne prone skin, occasional blemishes, comedones, combination skin (oily T-zone) | MODERATE ACNE | H | | H | H | 20 MINUTES | 144 |
| ACNE: MODERATE TO SEVERE | Persistent Acne and breakouts, red and inflamed pustules, facial redness, oily skin, visible pores | SEVERE ACNE | H | | | H | 20 MINUTES | 114 |
| REDNESS: VASCULAR | Visible blood vessels, environmental damage, thin and vulnerable skin | REDNESS | | | M | M | 20 MINUTES | 71.4 |
| REDNESS: SKIN TONE | Persistent facial redness, flushing, diffused redness (rosy cheeks), visible capillaries | INFLAMMATORY SKIN | | | | M | 20 MINUTES | 37.2 |
| SENSITIVE SKIN: INFLAMMATORY | Compromised skin, easily aggressed, blotchy redness, thin and vulnerable skin, visible capillaries | INFLAMMATORY SKIN | | | | M | 20 MINUTES | 37.2 |
| SENSITIVE SKIN: PROBLEM CONDITIONS | General sensitive and aggressed skin conditions, dry, irritated and reactive skin, bacterial conditions | SEVERE ACNE | H | | | H | 20 MINUTES | 114 |
| PSORIASIS | Red, flaky and scaly patches, irritated and itchy sores, inflammatory lesions | PSORIASIS | | | M | H | 20 MINUTES | 134.4 |
| WOUND HEALING | Reactive indications, skin trauma, open wounds, burns, sunburn, red scars | WOUND HEALING | | | M | H | 20 MINUTES | 134.4 |
| PAIN | Musculoskeletal pain | PAIN | | | H | H | 20 MINUTES | 192 |

PRE-SET PROTOCOL POWER LEVELS

NOTE: In the pre-set treatment programmes, power levels can be manually adjusted by selecting pre-set treatment - pre-set mode - custom. The user will then be presented with the option of LOW / Medium / High. Selected power levels will adjust all wavelength settings. To select individual wavelength settings, this can be done in the CUSTOM programme.

ACNE:** Mild to moderate, blue only for 10 minutes, then red and NIR start to deliver all three wavelengths for the remaining 10 minutes. *REDNESS:** Vascular. NIR only for 10 minutes, then red starts to deliver red and NIR for the remaining 10 minutes.