EBA 200

Ⓓ	Bedienungsanleitung	13
Ⓔ	Operating Manual	48
Ⓕ	Mode d'emploi	83
Ⓘ	Istruzioni per l'uso	118

Andreas Hettich GmbH & Co. KG
Föhrenstraße 12, D-78532 Tuttlingen / Germany
Phone +49 (0)7461 / 705-0
Fax +49 (0)7461 / 705-1125
info@hettichlab.com, service@hettichlab.com
www.hettichlab.com

Single Registration Number:

DE-MF-000010680

© 2023 by Andreas Hettich GmbH & Co. KG

All rights reserved. No part of this publication may be reproduced without the prior written permission of the copyright owner.

Änderungen vorbehalten! , Modifications reserved! , Sous réserve de modifications ! , Con riserva di modifiche!

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We hereby declare under our responsibility
that the designated device:

Type of device	Small centrifuges
Name	EBA 200
Basic UDI-DI	4050674010024QC
GMDN	36465
Classification	Medical Device, class IIa (Annex VIII, Chapter III, Rule 3)
according to	Regulation (EU) 2017/745 Annex IX
Involved Notified Body	mdc medical device certification GmbH; CE 0483 Kriegerstraße 6; 70191 Stuttgart; Germany

and its accessories, which are listed in the rela-
ted technical documentation and whose confor-
mity has been assessed together with the device,
complies with the relevant provisions of the
Regulation (EU) 2017/745 on medical devices.

Intended use

The **EBA 200** centrifuge is a medical device
according to the Regulation on Medical Devices
(EU) 2017/745.

The device is used to separate whole blood or
blood components of human origin into its
component parts for further processing for
therapeutic treatment.
Suitable collection systems are deployed by the
user for this purpose. Preparation of the blood for
its testing for suitability is included. The user can
set each of the variable physical parameters
within the limits set by the device.

The centrifuge may only be used by qualified
personnel in closed laboratories. The centrifuge
is only intended for the use referred to above.
Intended use also includes observing all
instructions in the Operating Manual and
compliance with the required inspection and
maintenance work.

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Any other use or use beyond this is considered improper. Andreas Hettich GmbH & Co. KG shall not be liable for any damage arising from this.

The device also complies to the applicable provisions of the following European directives, ordinances and standards

- 2006/42/EC "Directive on machinery"
- 2014/30/EU "EMC Directive"
- 2014/35/EU „Low Voltage Directive“
- 2011/65/EC "RoHS Directive"
(without involvement of a notified body)
- (EC) 1907/2006 „Regulation on REACH“
(without involvement of a notified body)

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This declaration of conformity is valid from 10.03.2023 until 25.08.2027

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We hereby declare under our responsibility that the designated device:

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Type of device **Small centrifuges**
Name **EBA 200**
Basic UDI-DI **4050674010024QC**
GMDN **36465**
Classification **Medical Device, class IIa
(Annex VIII, Chapter III, Rule 3)**
according to **Regulation (EU) 2017/745
Annex IX**
Involved
Notified Body **mdc medical device certification
GmbH; CE 0483
Kriegerstraße 6; 70191
Stuttgart; Germany**

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and its accessories, which are listed in the related technical documentation and whose conformity has been assessed together with the device, complies with the relevant provisions of the Regulation (EU) 2017/745 on medical devices.

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Intended use

The **EBA 200** centrifuge is a medical device according to the Regulation on Medical Devices (EU) 2017/745.

The device is used to separate whole blood or blood components of human origin into its component parts for further processing for therapeutic treatment.

Suitable collection systems are deployed by the user for this purpose. Preparation of the blood for its testing for suitability is included. The user can set each of the variable physical parameters within the limits set by the device.

The centrifuge may only be used by qualified personnel in closed laboratories. The centrifuge is only intended for the use referred to above. Intended use also includes observing all instructions in the Operating Manual and compliance with the required inspection and maintenance work.

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Any other use or use beyond this is considered improper. Andreas Hettich GmbH & Co. KG shall not be liable for any damage arising from this.

The device also complies to the applicable provisions of the following European directives, ordinances and standards

- 2006/42/EC "Directive on machinery"
- 2014/30/EU "EMC Directive"
- 2014/35/EU „Low Voltage Directive"
- 2011/65/EC "RoHS Directive"
(without involvement of a notified body)
- (EC) 1907/2006 „Regulation on REACH"
(without involvement of a notified body)

Applied relevant standards

DIN EN ISO 13485:2016; DIN EN ISO 14971:2020;
DIN EN 61010-1:2011-07; DIN EN 61010-2-020:2017;
DIN EN 61010-2-101:2002; DIN EN 61010-2-011:2017

Tuttlingen, 10.03.2023

Klaus-Günter Eberle
Gérant, Chief Executive Officer

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Type of device **Small centrifuges**
Name **EBA 200**
Basic UDI-DI **4050674010024QC**
GMDN **36465**
Classification **Medical Device, class IIa (Annex VIII, Chapter III, Rule 3)**
according to **Regulation (EU) 2017/745 Annex IX**
Involved Notified Body **mdc medical device certification GmbH; CE 0483 Kriegerstraße 6; 70191 Stuttgart; Germany**

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Intended use

The **EBA 200** centrifuge is a medical device according to the Regulation on Medical Devices (EU) 2017/745.
The device is used to separate whole blood or blood components of human origin into its component parts for further processing for therapeutic treatment.
Suitable collection systems are deployed by the user for this purpose. Preparation of the blood for its testing for suitability is included. The user can set each of the variable physical parameters within the limits set by the device.
The centrifuge may only be used by qualified personnel in closed laboratories. The centrifuge is only intended for the use referred to above. Intended use also includes observing all instructions in the Operating Manual and compliance with the required inspection and maintenance work.

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Any other use or use beyond this is considered improper. Andreas Hettich GmbH & Co. KG shall not be liable for any damage arising from this. The device also complies to the applicable provisions of the following European directives, ordinances and standards

- 2006/42/EC "Directive on machinery"
- 2014/30/EU "EMC Directive"
- 2014/35/EU "Low Voltage Directive"
- 2011/65/EC "RoHS Directive"
(without involvement of a notified body)
- (EC) 1907/2006 "Regulation on REACH"
(without involvement of a notified body)

Applied relevant standards

DIN EN ISO 13485:2016; DIN EN ISO 14971:2020;
DIN EN 61010-1:2011-07; DIN EN 61010-2-020:2017;
DIN EN 61010-2-101:2002; DIN EN 61010-2-011:2017

Tuttlingen, 10.03.2023

Klaus-Günter Eberle
Amministratore delegato, Chief Executive Officer

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- **USA:** QSR, 21CFR 820 "CFR Title 21 - Food and Drugs: TITLE 21- FOOD AND DRUGS, CHAPTER I - FOOD AND DRUG ADMINISTRATION DEPARTMENT OF HEALTH AND HUMAN SERVICES, SUBCHAPTER H - MEDICAL DEVICES, Part 820 QUALITY SYSTEM REGULATIONS"
- **Kanada:** CMDR, SOR/98-282 "Medical Devices Regulations"

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Standards and regulations which apply to this device

The device is a high-end technical product. It is subject to extensive testing and certification procedures according to the following standards and regulations in their respectively valid version:

Electrical and mechanical safety for design and final testing:

Standard series: IEC 61010 (conform to standards of DIN EN 61010)

- IEC 61010-1 "Safety requirements for electrical equipment for measurement, control, and laboratory use - Part 1: General requirements" (Pollution Degree 2, Excess-voltage category II)
- IEC 61010-2-010 "Safety requirements for electrical equipment for measurement, control and laboratory use - Part 2-010: Particular requirements for laboratory equipment for the heating of materials" (applied to heated centrifuges only)
- IEC 61010-2-011 "Safety requirements for electrical equipment for measurement, control and laboratory use - Part 2-011: Particular requirements for refrigerating equipment" (only valid for centrifuges with cooling)
- IEC 61010-2-020 "Safety requirements for electrical equipment for measurement, control, and laboratory use - Part 2-020: Particular requirements for laboratory centrifuges"
- IEC 61010-2-101 "Safety requirements for electrical equipment for measurement, control and laboratory use - Part 2-101: Particular requirements for in vitro diagnostic (IVD) medical equipment"

Electromagnetic Compatibility:

- EN 61326-1 "Electrical equipment for measurement, control and laboratory use - EMC requirements - Part 1: General requirements"

Risk management:

- DIN EN ISO 14971 "Application of risk management to medical devices"

Restriction of Hazardous Substances (RoHS II):

- EN 50581 "Technical documentation for assessing electric and electronic devices with regard to the restriction of hazardous substances"

European Directives applicable to conformity assessment procedures:

- EU Quality Management Certificate in accordance with Regulation (EU) 2017/745 on medical devices Annex IX, Chapter I (Quality Management System)
The certificate is issued by the following notified body:

mdc medical device certification GmbH – Notified Body CE 0483

Tel: +49 (0)711 253597 0

Fax: +49 (0)711 253597 10

Email: mdc@mdc-ce.de

Website: www.mdc-ce.de

Post: Kriegerstrasse 6, D-70191 Stuttgart, Germany

Monitoring is carried out in accordance with Annex IX, Section 3 of Regulation (EU) 2017/745.

Risk Class IIa medical device according to Annex VIII, Rule 3 of Regulation (EU) 2017/745.

- Directive 2011/65/EU for the restriction of use of certain hazardous substances in electric and electronic devices. Carrying out the EC conformity assessment process is the sole responsibility of the manufacturer, without the involvement of a notified body.

Applied medical device regulations outside Europe:

- **USA:** QSR, 21CFR 820 "CFR Title 21 - Food and Drugs: TITLE 21- FOOD AND DRUGS, CHAPTER I - FOOD AND DRUG ADMINISTRATION DEPARTMENT OF HEALTH AND HUMAN SERVICES, SUBCHAPTER H - MEDICAL DEVICES, Part 820 QUALITY SYSTEM REGULATIONS"
- **Canada:** CMDR, SOR/98-282 "Medical Devices Regulations"

Certified quality management system according to

- ISO 9001 "Quality management systems – Requirements"
- ISO13485 "Medical devices - Quality management systems - Requirements for regulatory purposes"

Environmental management system according to

- ISO 14001 "Environmental management systems - Requirements with guidance for use"

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