
- de** Gebrauchsanweisung
- en** User's Manual
- fr** Mode d'emploi
- it** Istruzioni per l'uso
- es** Instrucciones de empleo
- pt** Manual de operação
- ru** Руководство по применению

Stimmgabeln


QM-1-027F

de	Stimmungsgabeln	4
en	Tuning Forks	8
fr	Diapasons	12
it	Diapason	16
es	Diapasones	20
pt	Diapasões	24
ru	Камертоны	28

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
User's Manual Tuning Forks

e.g. Rydel Seiffer-type tuning fork,
c 128/C 64 Hz with damper
and plastic foot

adjustable

e.g. Hartmann-type
tuning fork, c2 512 Hz
without damper and plastic foot

e.g. Lucae-type tuning fork,
c2 512 Hz without damper,
with plastic foot

Dear Customer, thank you for choosing a KaWe product. Our products are known for their high quality and long life. This KaWe product meets EC Standards 93/42/EEC (standards for medical products).

Please read these instructions thoroughly and carefully before attempting to use this product and heed the given care instructions.

Familiarize yourself fully with this product before using it.

Use: These tuning forks are only to be used by authorized, trained personnel. The tuning forks create sound that is then conducted through the air and bone. The tuning forks are only to be placed on healthy skin. The plastic foot or end of the vibrating fork is placed on to the part of the body that is being examined. These are primarily spots on the extremities where the joints are not covered by muscle or fatty tissue. In order to ensure long-term monitoring, the results of each inspected area should each be compared with one another over time. For this purpose, the response of the patient is needed in order to determine at which point of the fading process of the fork's oscillation the patient no longer feels the vibration.

Use of the Rydel Seiffer-type tuning forks: These tuning forks are equipped with 2 scaled dampers. The scale consists of a contoured and an out-lined arrow with 8 numbered, lateral division lines. When using, pay attention to the adjustment of the lower edge of the damper on the marked lines of both prongs. The fork is held in one hand at the stem and is vibrated by striking the adjusting screw against the palm of the other hand. The tuning fork can also be vibrated by plucking the prongs with the fingers. When vibrated, each arrow appears as two triangles, the sides of which intersect at a point that travels upward as the vibration amplitude decreases. The plastic foot or end of the vibrating fork is then placed on to the part of the body that is being examined. These are primarily spots on the extremities where the joints are not covered by muscle or fatty tissue. In order to ensure long-term monitoring, the results of each inspected area should each be compared with each other over time. For this purpose, the response of the patient is needed in order to determine at which point of the fading process of the fork's oscillation the patient no longer feels the vibration. At this point the reading can be taken using the lateral markings. This value is to be compared with previous measurements. Fluctuations in readings for individual patients are possible. Should for example, a patient sense vibrations up to the 6th division line, the value to be recorded is 6/8.

As is known, one can detect a beginning deterioration of the nervous system through a decrease in vibration perception. A replication of the magnitude of these vibrations is a prerequisite for the practical application of this diagnosis.

Purpose: When hearing, normally airborne sound travels through the outer ear channel to the eardrum. The signal is transferred by the ossicle of the middle ear to the inner ear and is there taken in by the receptor cells (aerotympanal conduction). Sound oscillations can also be transferred through the bones, by-passing the middle ear, directly to the inner ear (bone conduction). The transfer of sound oscillations through air channels or bone channels depends on the frequency and differs between aerotympanal conduction and bone conduction.

Unsuited use/contraindication of the product: Any use other than that described here is not in accordance with the intended use of the unit. The manufacturer is not liable for any resulting damages. The user alone bears the risk.

Further information, maintenance, storage: When used and stored properly, this product will serve you many years.

Guarantee: When used under normal circumstances and with attention to our instructions (the User's Manual) we guarantee this product for two years after the date of purchase. Should you need further information or should your instrument require repair, please contact your dealer!

Warning: Please keep in mind that even minor damage to the tuning fork can influence the oscillation frequency.

Limitation of reuse preparation: Frequent re-preparation has only a minimal effect on the instrument. The end of the product's life is normally due to wear and tear and damage by use. It is recommended that these instruments be re-prepared as soon as possible after their use.

Instructions

Storage and transport: No special requirements.

Cleaning preparation: No special requirements.

Manual cleaning: Clean the tuning forks by wiping them with a damp cloth or surface cleaning agent. Heed hereby the concentrations given by the manufacturer.

Automatic cleaning: Please ensure that the cleaning system is suitable for medical instruments.

Only cleaning agents suitable for metals are to be used. The concentrations and reaction times given by the manufacturer of the cleaning agents are to be followed.

Disinfection: Parts that come into contact with the patient are to be disinfected after each use. Alcohol-based surface disinfecting agents are to be used. Please heed hereby the instruc-

tions and reaction times given by the manufacturer. Especially the surfaces of the instruments that touch the patient are to be re-prepared each time the instrument is used.

Maintenance: Damaged parts must be removed and/or replaced.

Check and performance test: Check moving parts such as the dampers for smooth operation. Check for damage and wear. Before use, ensure that all screws are tightened.

Packaging: packaged individually: a standard polyethylene bag can be used. The size of the bag must correspond to the size of the instrument so that it closes properly without straining the bag closure.

Sterilization: No sterilization required, as it is not a surgical instrument.

Storage: No special requirements.

Disposal: The product can be disposed of at a location that accepts scrap metal.

Common accessories: Further information about these articles can be found on our homepage: www.kawemed.de.

Manufacturer: **KaWe**


Contacting the manufacturer: Write to or phone your dealer or dial +49-7141-68188-0.

Symbol key:

	Manufacturer		Caution!
	Date of manufacturer		Complies with relevant EU guidelines
	Heed the User's Manual		

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KIRCHNER & WILHELM GmbH + Co. KG