









# Prontosan®

## BREAKS THE BIOFILM CYCLE

A proactive approach using a combination strategy of Prontosan Solution and Prontosan Gel X as part of wound bed preparation may prove helpful and aims to:

- ☐ Reduce the biofilm burden (**Prontosan Solution**)
- ☐ Prevent reconstitution of the biofilm (**Prontosan Gel X**)

### Appropriate Time for Appropriate Wound

DESCRIPTION OF WOUND	OBJECTIVE	HOW TO USE
ACUTE WOUND - SURGICAL PRIMARY & SECONDARY INTENTION HEALING		
<input type="checkbox"/> High risk patient* <input type="checkbox"/> No slough <input type="checkbox"/> Minimal exudate	<input type="checkbox"/> Cleans <input type="checkbox"/> Prevents biofilm / complications	Rinse with Solution 
ACUTE WOUND e.g. trauma		
<input type="checkbox"/> Debris <input type="checkbox"/> Haematoma	<input type="checkbox"/> Cleans <input type="checkbox"/> Prevents biofilm / complications	Soak with Solution 
CHRONIC WOUND – GRANULATING		
<input type="checkbox"/> High risk patient* <input type="checkbox"/> Low exudate	<input type="checkbox"/> Cleans <input type="checkbox"/> Prevents biofilm / complications	Soak with Solution Consider Gel X  + 
CHRONIC WOUND		
<input type="checkbox"/> Light slough <input type="checkbox"/> Low exudate	<input type="checkbox"/> Cleans <input type="checkbox"/> Prevents biofilm / complications	Soak with Solution Apply Gel X  + 
CHRONIC WOUND - CRITICALLY COLONISED/INFECTED		
<input type="checkbox"/> Medium / high exudate <input type="checkbox"/> Static wound <input type="checkbox"/> Slough	<input type="checkbox"/> Cleans <input type="checkbox"/> Prevents biofilm / complications	Soak with Solution Apply Gel X  + 

\*High risk patient: Co-morbidities such as Diabetes, immuno-compromised, steroidal use, patients with previous wound infection and or biofilm and slough.

Compatible with other dressings | Once opened, Prontosan Solution and Gel X can be used for up to 8 weeks (single patient use)