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
² Steris Device Testing (2016) "External Summary Report Device Testing" Dated 20 October 2016. Report for Vitalograph Ireland Ltd. Report Number: 10250035. A Mentor, OH. STERIS Corporation.

³ Vitalograph (Ireland) Ltd (2014) "Testing the recommended cleaning method – Report Pneumotrac Flowhead" Dated 13 May 2014. Internal Vitalograph report. Report Number: P229_0038.

⁴ Vitalograph (Ireland) Ltd (2015) "Testing the recommended cleaning method using PeraSafe - Report Pneumotrac Flowhead" Dated 20 Apr 2015. Internal Vitalograph report. Report Number: P229_0040

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60062

Referenzen

1. Steris Device Testing (2016) "External Summary Report Device Testing" Dated 20 October 2016. Internal report for Vitalograph Ireland Ltd. Report Number: 10250035. A Mentor, OH. STERIS Corporation.*
2. Vitalograph (Ireland) Ltd (2014) "Testing the recommended cleaning method – Report Pneumotrac Flowhead" Dated 13 May 2014. Internal Vitalograph report. Report number: P229_0038.*

3. Vitalograph (2019), "Hygiene Policy". Internal Vitalograph policy. Document number: SOP_0523.*

Literaturverzeichnis

1. 1. Bentz, J. R. (2019). "Bacterial Filtration Efficiency (BFE) at an Increased Challenge Level Final" Study Number 1138681-SGI; "Viral Filtration Efficiency (VFE) at an Increased Challenge Level Final" Study Number 1138680-S01. Internal reports for Vitalograph Ireland Ltd dated 17 January 2019. Nelson Laboratories, Salt Lake City.*
2. Bracci, M. et al (2011). "Risk of bacterial cross infection associated with inspiration through flowbased spirometers." American Journal of Infection Control 39(1): 50-55. DOI: <https://doi.org/10.1016/j.ajic.2010.04.215>.
3. Dunne, C (2019). "Calibrated Flow Bioburden testing of Vitalograph Alpha Flow Heads" Dated 16 July 2019. Internal report for Vitalograph Ireland Ltd.*
4. FDA (2015) "Reprocessing Medical Devices in Health Care Settings: Validation Methods and Labeling Guidance for Industry and Food and Drug Administration Staff" U.S. Food & Drug Administration.
5. FDA (2018). "What are Reusable Medical Devices?" U.S. Food & Drug Administration. Accessed 4 July 2019, from <https://www.fda.gov/medical-devices/reprocessing-reusable-medical-devices/what-are-reusable-medical-devices>
6. Kendrick, A. H. et al (2003). "Infection control of lung function equipment: a practical approach." Respiratory Medicine 97(11): 1163-1179. DOI: [https://doi.org/10.1016/S0954-6111\(03\)00223-3](https://doi.org/10.1016/S0954-6111(03)00223-3)
7. Loveday, H. P. et al (2014). "epic3: National Evidence-Based Guidelines for Preventing Healthcare-Associated Infections in NHS Hospitals in England." Journal of Hospital Infection 86: S1-S70. DOI: [https://doi.org/10.1016/S0195-6701\(13\)60012-2](https://doi.org/10.1016/S0195-6701(13)60012-2).
8. NHS (2017). Community Infection Prevention and Control Guidance for General Practice. Infection Prevention Control. UK. Accessed 28 May 2019, from <https://www.infectionpreventioncontrol.co.uk/>
9. Rutala, W. A., D. J. Weber and HICPAC (2008). "Guideline for Disinfection and Sterilization in Healthcare Facilities." CDC Infection Control Accessed 28 May 2019, from <https://www.cdc.gov/infectioncontrol/guidelines>
*Daten aus internen Berichten und Richtlinien können durch Vitalograph auf Anfrage zur Verfügung gestellt werden.