bosotherm diagnostic

- (D)** **Kontaktloses IR-Thermometer**
Gebrauchsanweisung
- (GB)** **Non-Contact IR Thermometer**
User Instructions
- (F)** **Thermomètre infrarouge sans contact**
Mode d'emploi
- (I)** **Termometro IR senza contatto**
Istruzioni per l'uso
- (E)** **Termómetro IR exento de contacto**
Manual del usuario

Kontaktloses IR-Thermometer bosotherm diagnostic

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Non-Contact IR Thermometer bosotherm diagnostic

User Instructions

Congratulations on your purchase of the boso infrared thermometer, bosotherm diagnostic. This digital, infrared thermometer is a quality, state of the art product.

With its unique technology, the bosotherm diagnostic enables stable measurement values with every measurement which are free of any disruptive thermal influences. The device performs an auto-test each time it is activated in order to ensure the specified measurement accuracy.

Intended Use

Determination of the body temperature of people through non-contact measurement of the surface temperature of the forehead, as well as the determination of the temperature of objects.

Carefully read these instructions before using the device and store them in a secure location.

In order to get help with commissioning, use or maintenance, please contact your retailer or the manufacturer.

An unexpected operating state or an occurrence which worsens or could worsen a health condition must be reported to the manufacturer immediately.

By affixing the CE marking, the manufacturer confirms that the product is in conformity with the applicable European Medical Device Directive 93/42/EEC. The four-digit identification number identifies the notified body.


Contents

Device Overview / Display	26
Scope of Delivery	27
Important Safety Instructions	28
Device Properties	31
Operation of the Device	32
Error Messages and Troubleshooting	38
Cleaning and Maintenance	40
Warranty / Customer Service	42
Technical Data	43
Disposal	45


The device must be installed and commissioned according to the information in these user instructions.

Wireless communication devices, such as wireless home networking devices, cell phones, cordless telephones and their base stations and walkie-talkies can influence this thermometer. Therefore, a clearance of at least 3.3 meters should be maintained from such devices.

Device Overview

Display

Scope of Supply

- 1 IR Thermometer
bosotherm diagnostic

- 1 Storage bag

- 2 LR03 batteries (AAA)
(inserted in device)

- 1 User instructions

Important Safety Instructions

Attention

- This device cannot replace a doctor's visit or a medical diagnosis.
- Please be aware of the risk of burns after heating objects in the microwave, since a surface measurement can produce a lower temperature than the temperature inside the object.
- Do not use the device if you believe that it is damaged or if you have noticed something unusual.
- The device is not waterproof – do not immerse it in water or other liquids.
- Ensure that this device is not accessible to children in order to prevent swallowing of batteries or the battery compartment cover.
- This device may not be opened, dismantled or modified.
- Before use, the device should be acclimated to the ambient temperature (that is, room temperature). If a change in ambient temperature should occur, such as due to a temperature difference between the storage location and the usage location of the device, you should wait at least 30 minutes before use so that the device can adapt to a stable ambient temperature (usage location)!
- Soiling on the IR sensor must be carefully removed with a dry cloth or a cotton swab. Do not use paper towels or paper tissues to clean the sensor.
- The device contains sensitive electronic components – in order to ensure its capability to perform, it must be used and stored according to the descriptions in the user instructions.

Protect the device from

- Extreme temperatures
 - Shocks and vibrations
 - Dust
 - Direct sunlight
 - Extreme humidity
- Please remove the batteries if the device is not used for an extended period of time.
 - This device measures the skin temperature on the forehead in order to determine human body temperature. Skin temperature can vary in dependence on various influences, such as due to bathing, swimming or other athletic activities, due to sweating, due to the use of ice packs and the extended effects of high or low temperatures (outdoors). Before taking temperature measurements, please wait at least 30 minutes so that the body can adapt to a stable ambient temperature (room temperature).
 - Water, other surface liquids (e.g. sweat) or condensate can negatively influence infrared measurements.
 - Please ensure that the skin surface or the surface of the object which is measured, as well as the sensor lens, are dry before performing a temperature measurement, in order to ensure the accuracy of this measurement.
 - In the initial stage of a febrile illness, a physiological phenomenon called vasoconstriction (vascular constriction) can occur, which leads to a lower skin temperature. For this reason, a measurement with this thermometer may return an unusually low result.
 - If the measurement result should not correspond with the condition of the patient or should be unusually low, the measurement should be repeated at 15 minute intervals or compared to a different measurement of core body temperature.
 - Do not use the thermometer immediately after the following:
 - Drinking or Eating
 - Exercising
 - Nursing a baby

Important Notes about Taking Temperature Measurement of Infants & Juveniles!

Doctors recommend rectal measurement for newborn infants within the first 6 months, as all other measuring methods might lead to ambiguous results. If using a non-contact thermometer on infants of this age range, we always recommend verifying the temperature readings with a rectal measurement!

In the following situation it is recommended that three temperatures are taken with the highest one taken as the reading:

1. New born infants in the first 100 days.
2. Children under age of three with a compromised immune system and for whom the presence or absence of fever is critical.
3. When the user is learning how to use the thermometer for the first time, until he/she has familiarized himself/herself with the device and capable of obtaining consistent readings.
4. If the measurement is surprisingly low.

Note

Readings from different measurement sites should not be compared because the normal body temperature varies between different measurement sites and the time of the day, being highest in the evening and lowest about one hour before waking up.

Normal body temperature ranges:

- Axillar: 34.7 - 37.3°C
- Oral: 35.5 - 37.5°C
- Rectal: 36.6 - 38.0°C
- bosotherm diagnostic: 35.4 - 37.4°C
- Objects with a low IR emission level can return temperature data which is below the actual object temperature.
- Do not touch the lens of the infrared sensor – contaminations on the lens can cause inaccuracies in the temperature display.

Device Properties

Temperature measurements with the bosotherm diagnostic

The bosotherm diagnostic thermometer uses a sensitive sensor technology in order to receive the emission of infrared radiation from the human body or objects and, based on this, to perform precise measurements of the human body temperature or the temperature of objects.

Safe and comfortable

- No direct contact! Temperature measurement without skin contact
- No danger of glass breakage or swallowing mercury

Fast measurement

- Measurements within 1 second

Various measurement methods

- The device enables 2 different methods for measuring 1. human body temperature and 2. the temperature of objects

Saving and recall of temperature measurement

- 30 temperature measurements can be saved and recalled with the device.
- Direct access to the recorded temperature measurements of the body or of an object.

Fever alarm*

- At temperature measurements of $\geq 37.5^{\circ}\text{C}$, an acoustic signal automatically sounds in order to indicate an elevated body temperature.

* only for measurements on the human body

Operation of the Device

Before using for the first time

The bosotherm diagnostic device contains 2 X LR03 (AAA) alkaline batteries which are equipped with an isolating strip. This isolating strip must be removed from the battery compartment before using the bosotherm diagnostic for the first time.

Measurement of Human Body Temperature

1. Press the  button in order to start the bosotherm diagnostic in body temperature mode. After activating, the device performs a self-test and switches to body temperature mode after a few seconds. When the bosotherm diagnostic is ready for a measurement in body temperature mode, "- - -" blinks in the display and the symbol for body temperature mode appears.

2. To measure the temperature of the human body, hold the device at a distance of about 5 cm (or less) in front of the center of the forehead so that the sensor points toward the forehead.

3. Once the bosotherm diagnostic is in the correct position, press the START button to perform a measurement of the body temperature. The measurement process takes about one second. The display indicates that the measurement is running.

4. If the measured body temperature is below 37.5°C , the measurement result is displayed with a long signal tone.

5. If the measured body temperature is at or over 37.5°C (as is the case with a fever), there are 3 short signal tones.

6. After performing a temperature measurement, the bosotherm diagnostic shows the measurement result continuously, and switches off after 1 minute. As soon as the device is ready for a new measurement, there is a short signal tone and the " $^{\circ}\text{C}$ " symbol begins to blink. To perform another temperature measurement, press the START button and the bosotherm diagnostic starts a new measurement. Once the new measurement is complete, there is a short signal tone and the new measurement result is displayed.

→ Push START button

7. If the measurement returns a result that is above or below the specific temperature range for this mode, "H" or "L" is shown in the display and three short signal tones are emitted. After about 5 seconds, the bosotherm diagnostic is ready to measure again.

or

→

8. Press the  button in order to switch off the bosotherm diagnostic. The device automatically switches off after 1 minute.

→ Push  button

→


Recall of Last Temperature Measurement In Body Temperature Mode

To recall the latest temperature measurement in body temperature mode, press the M button while the bosotherm diagnostic is in body temperature mode. The last body temperature measurement is displayed together with the blinking M symbol for about 5 seconds.

Temperature Measurements of Objects

1. Press the START button and hold it while activating the bosotherm diagnostic with the  button. After activating, the device performs a self-test and will be in object mode after a few seconds. When the bosotherm diagnostic is ready for a temperature measurement of an object, the " - - - " symbol blinks in the display and the symbol for object mode appears.

2. To correctly measure the temperature of an object, hold the device at a distance of about 5 cm (or less) in front of the surface so that the sensor points toward the object.

3. Once the bosotherm diagnostic is in the correct position, press the START button to perform a measurement of the object's temperature. The measurement process takes about one second. The display indicates that the measurement is running.

4. Once the measurement is complete, there is a long signal tone and the measurement result is displayed.

5. After performing a temperature measurement, the bosotherm diagnostic shows the measurement result continuously, and switches off after 1 minute. As soon as the device is ready for a new measurement, there is a short signal tone and the "°C" symbol begins to blink. To perform another temperature measurement, press the START button and the bosotherm diagnostic starts a new measurement.

→ Push START button

6. If the measurement returns a result that is above or below the specific temperature range for this mode, "H" or "L" is shown in the display and three short signal tones are emitted. After about 5 seconds, the bosotherm diagnostic is ready to measure again.

or

→

7. Press the  button in order to switch off the bosotherm diagnostic. The device automatically switches off after 1 minute.

→ Push  button

→


Recall of Last Temperature Measurement In Object Mode

To recall the latest temperature measurement in object mode, press the M button while the bosotherm diagnostic is in object mode. The last temperature measurement is displayed together with the blinking M symbol for about 5 seconds.


Recall All Temperature Measurements In Memory Mode

1. Press the M button while device is deactivated in order to activate the memory mode. If no measurement results are saved, the bosotherm diagnostic displays "0." If measurement results are saved, the display begins with the latest temperature measurement. "1" appears in the display and automatically changes to display the temperature measurement.

2. Press the M button multiple times in succession to switch to a specific recorded measurement result.

3. Press the  button in order to switch off the device. The bosotherm diagnostic switches off automatically after 30 seconds.

Clear Measurement Value Memory

To delete all temperature measurements in the memory mode of the bosotherm diagnostic, hold the START button in memory mode for 5 seconds or longer. The bosotherm diagnostic displays "CL" and all saved temperature measurements are deleted. The display then changes to "0."

Error Messages and Troubleshooting

Self-test

- If the bosotherm diagnostic determines an error during the self-test phase after activation, "Er" is displayed and 3 short acoustic signals are emitted; the device then switches off after 10 seconds.

Troubleshooting:

Remove the batteries from the device and wait a few minutes before replacing them in the device.

Attempt to activate the device again. If the error message should appear again, contact a local dealer.

Low batteries display

- If the batteries are low, the battery symbol appears. Please change the batteries as quickly as possible.

- If the voltage falls below the value required for the function, the bosotherm diagnostic solely displays a battery symbol after pressing the  button.

- If the batteries are completely discharged, no display appears after pressing the  button.

Troubleshooting:

Replace the batteries!

Measurement outside of specified ambient conditions

- The bosotherm diagnostic automatically measures the ambient temperature in body temperature mode and object mode. If the ambient temperature is below or above the defined temperature range, a warning is displayed and the temperature measurement is not carried out.
- If the ambient temperature (e.g. room temperature) is $>40^{\circ}\text{C}$, the warning notice "Er.H" appears and there are 3 short acoustic signals. The bosotherm diagnostic switches off automatically after 10 seconds.

- If the ambient temperature (e.g. room temperature) is low in body temperature mode (ambient temperature $<15^{\circ}\text{C}$) or in object mode (ambient temperature $<5^{\circ}\text{C}$), the warning notice "Er.L" appears and there are 3 short acoustic signals. The bosotherm diagnostic switches off automatically after 10 seconds.

Troubleshooting:

Please only use the bosotherm diagnostic under the specified environmental conditions. Move the device to an environment with the specified conditions and wait 30 minutes before using it.

Cleaning and Maintenance

Storage of the bosotherm diagnostic

- The bosotherm diagnostic contains highly sensitive electronic components. In order to ensure appropriate performance and long service life, the bosotherm diagnostic should be stored in a dry location, out of direct sunlight, and protected from shocks and vibrations.
- Avoid storing in a wet location, or in the vicinity of heat sources.
- Please remove the batteries from the bosotherm diagnostic before storing for an extended period.

Cleaning and disinfection of the bosotherm diagnostic

- Use an alcohol-soaked swab or a cotton cloth moistened with alcohol (70% isopropyl) for cleaning. For wiping disinfection (exposure time min. 5 minutes) of the device, we recommend the disinfectant Microzid Liquid Sensitiv (Schülke & Mayr).
- Allow the bosotherm diagnostic to dry completely before using again and before performing more temperature measurements.
- The device is not waterproof – do not immerse it in water or other cleaning fluids.
- Do not use aggressive or corrosive cleaning agents or solutions - this can damage the device.

Changing batteries

- The bosotherm diagnostic is operated by two LR03 (AAA) alkaline batteries. Please replace the old batteries with new ones when the battery symbol blinks.
- Open the battery compartment cover and remove the batteries from the battery compartment. Insert new LR03 batteries and take note of the symbols for the poles in the battery compartment: The (+) pole must correspond with the "+" symbol and the (-) pole with the "-" symbol. Close the battery compartment again after inserting the batteries.
- Saved temperature measurements are retained during a battery change.
- Avoid using new and old batteries from different manufacturers.
- Do not leave old batteries in the bosotherm diagnostic device. This can cause the batteries to leak and can damage the device.

Regular Checks on the Accuracy

Regular checks on the accuracy of the measurement system should be carried out in accordance with the valid directive of the individual countries, governing this issue.

Warranty / Customer Service

We provide a factory warranty for 2 years from the date of purchase for this product. The purchase date must be substantiated by receipt.

Within the warranty period, defects resulting from material or production errors will be resolved free of charge.

As a result of warranty service, no extension of the warranty period is applied to the entire device, but rather only to the replaced components.

Abrasion due to wear, transport damages as well as any damages caused by improper use (e.g. failure to comply with usage instructions) or by the actions of unauthorized parties are excluded from the warranty coverage. No damage compensation claims against us are justified by the warranty.

Legal defect claims of the purchaser according to § 437 BGB are not restricted. In case of warranty claim, send the device with the original purchase receipt to:

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Technical Data

Measuring range:

- Range for measurement of body temperature 34°C ~ 42.2°C
- Range for measurement of the temperature of objects: 0.1°C ~ 99.9°C

Resolution: 0.1°C

Measurement accuracy under laboratory reference conditions:

- Accuracy when measuring body temperature:
±0.2°C in the range from 35°C ~ 42°C
±0.3°C in the range < 35°C and > 42°C
- Accuracy when measuring the temperature of objects: ±1.0°C in the range from 0.1°C ~ 99.9°C

Measurement distance: max. 5 cm

Display LCD

Memory: 30 measured values

Backlight: blue

Operating conditions:

- Temperature range: 15°C ~ 40°C
- Humidity range: up to 95% relative humidity (no condensation)

Storage and transport:

- Temperature range: -20°C ~ +55°C
- Humidity range: up to 95% relative humidity (no condensation)

Automatic shut-off:

- 60 seconds after completing a temperature measurement
- 30 seconds in memory mode
- 10 seconds after error message and in case of deviation in ambient temperature from defined temperature range

Voltage/batteries: 1.5V x 2 AAA (LR03)

Applied part: Type BF ()

Service life of batteries: at least 1,000 measurements

Anticipated service life of the device: 10 Years

Dimensions: approx. 140 mm (L) x 43 mm (W) x 47 mm (H)

Weight: approx. 93 g incl. batteries

Protection against penetration by water or solid materials: IP 22

Reference Standards: ASTM E 1965-98
IEC 60601-1
IEC 60601-1-2
IEC 60601-1-11
ISO 80601-2-56

The bosotherm diagnostic conforms to the provisions of the EU Guidelines 93/42/EEC for Class IIa medical devices.

Technical changes reserved.

Disposal

Batteries and electronic devices must be disposed of in accordance with local regulations. Please do not discard batteries or electronic devices in household trash.

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BOSCH + SOHN GmbH u. Co. KG