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Englisch

Information on the variousmodels of Riester aneroid sphgmomanometers

You have acquired a high-quality **Riester** precision aneroid sphgmomanometer, which has been manufactured in accordance with the ISO 81060-1: 2007 standard [DIN EN ISO 81060-1: 2013] and is subject to the strictest quality controls. The excellent quality will guarantee years of reliable measurements.

Intended purpose: The aneroid sphgmomanometers from **Riester** are used by physicians and persons trained in auscultatory blood pressure measurement to determine the systolic and diastolic blood pressure in human beings and neonates. **Riester** aneroid sphgmomanometers are not intended exclusively for blood pressure measurement on healthy skin on the upper arm or thigh. The product is normally used professionally in a medical practice or a hospital. The sphgmomanometer/aneroid manometer is intended for use as a diagnostic aid.

Product description / indication: All blood pressure measuring devices / aneroid manometers distributed by company **Rudolf Riester GmbH** are of the same basic structure. Blood pressure measuring devices / aneroid manometers consist of a measuring system, cuff and inflation system (consisting of the pump (ball with non-return valve), and the screw valve (air lock or air drain)). All devices work according to the same measuring principle: Indirect blood pressure measurement according to the Korotkoff method.

Warnings / Contraindications:

–If there is overpressure in the cuff, you can reduce the pressure with the quick release.

For units with air release valve: fully open the air release screw.

For devices with push-button valve: Press push-button valve completely.

Use the ball to pump up the cuff up approximately 20mmHg above the expected systolic blood pressure value (= the upper value). Never inflate higher than 300mmHg.

The cuffs are offered with latex-free and latex materials. These are indicated by corresponding symbols on the cuff.

–Pay attention to the paragraph „blood pressure measurement“ in the instructions for use, the blood pressure measurement can only done from doctors and persons trained in auscultatory blood pressure measurement.

Observe the measurement and storage conditions in the instructions for use under „Technical data“ and the information on the label.

–Check the pointer zero position as described in „Test of precision“.

For Germany, the Medical Device Ordinance applies. See also point „metrological control“.

–A description of where the LOT No. is located can be found in the cuff table at the end of the instructions for use.

1. Models exacta® and sphgmotensiophone
They are equipped with 2-tube cuffs because the pressure gauge is not directly connected to the inflation system.

2. Models e-mega, R1 shock-proof®, minimus® II, precisa® N, precisa® N shock-proof, babyphon® and **ri-sani® 1-tube.**
They are equipped with 1-tube cuffs. Here the inflation system is directly connected to the manometer.

3. Models e-mega, R1 shock-proof®, minimus® II, precisa® N, precisa® N shock-proof, babyphon® and **ri-sani® 1-tube.**
They are equipped with double tube cuffs. Here the inflation system is directly connected to the manometer.

4. Models sanaphon® and ri-sani®
These devices are devices for self-measurement. The manometer and the inflation system form a unit. In the case of the ri-sani® model, the chestpiece is integrated, which absorbs the Korotkoff sounds and forwards them to the ear via the enclosed stethoscope, which must be screwed into the thread on the outside of the cuff.

5. Models bigben Round / Square desk, wall, mobile and anesthesia and ri-former® bigben
These cuffs with push-button valve. Press push-button valve completely. Use the ball to pump up the cuff up approximately 20mmHg above the expected systolic blood pressure value (= the upper value). Never inflate higher than 300mmHg.

The cuffs are offered with latex-free and latex materials. These are indicated by corresponding symbols on the cuff.

Installation instructions for the different models (with the exception of the desk model):

5.1. Wall model / ri-former® bigben
Remove the wall nut underneath the cuff basket and take off the wall bracket. Hold the bracket against the wall at the desired location and mark the positions of the mounting holes for drilling. Hold the holes and insert the wall bracket in the wall and tighten the screws. Place the unit on the wall bracket so that the top part of the wall bracket engages the edge of the cuff basket and the bottom part fits over the screw that protrudes from the underside of the cuff basket. Now replace and tighten the wing nut on the protruding screw.

5.2. mobile stand version / bigben

5.3. mobile stand version / bigben
1. Please assemble the mobile stand first and thereby take note of the installation instruction attached.
2. After assembling the mobile stand, please screw the device (clockwise) into the mobile stand.

Adjustment:
By opening the retaining screw the desired height can be adjusted. After adjusting the height the mobile stand must be fastened by the retaining screw again.

5.3. anesthesia model / bigben

Remove the wing nut underneath the cuff basket and take off the wall bracket. Affix the wall bracket to the back of 100% rigid plastic bracket so that the wall bracket is attached to the wall. Place the unit on the wall bracket so that the top part of the wall bracket engages the edge of the cuff basket and the bottom part fits over the screw that protrudes from the underside of the cuff basket. Now replace and tighten the wing nut on the protruding screw.

again filling up with blood.

- The blood pressure measurement has been completed.
- We would like to point out to you that a unit intended for home use does not replace regular visits to the doctor and that only the doctor can accurately analyze your measured values.

How to care for the aneroid sphgmomanometer

General information
The goal of cleaning and disinfection of medical products is the protection of patients, users and third persons and conserving the value of the medical products. On account of the product design and the used material, no fixed limit of maximum permissible cleaning cycles can be fixed. The lifetime of the medical products depends on their function and on an appropriate treatment of the devices. Before returning faulty products for repair they must have gone through the described reprocessing process.

1. Manometer and bulb
Manometer and bulb can be cleaned outside with a humid cloth until optical cleanliness is given.

ATTENTION!
Never place the manometer in liquid! This item is not approved for automated reprocessing and sterilization. These procedures cause irreparable damage!

2. Cuffs
Cotton and Nylon velcro cuff (latex and latex free)
Cleaning:
After removing the bladder, wipe the nylon-velcro covers with a damp cloth. Alternatively, these can be washed with soap and cold water like all the other cuffs. If you decide on the latter course, rinse the cuffs with clear water afterwards and let them air dry. Wipe the bladder and tubes with a damp cloth.

Disinfection:
After removing the bladder, wash the cuff covers in cold water to which disinfectant has been added. Afterwards, the cuffs must be dried only with clean, moistened cloths. The lifetime of the medical products depends on their function and on a appropriate treatment of the devices. Before returning faulty products for repair they must have gone through the described reprocessing process.

Disinfectable one piece cuff
Cleaning:
The cuff can be wiped with a damp cloth. Alternatively, it can be washed with soap and cold water like all the other cuffs. Please rinse the cuff with clear water afterwards. In addition this cuff can be washed up to 60° C in the washing machine, please ensure that no liquid is remaining in the cuff. This can affect the measurement results negatively and it can damage the manometer technology.

Disinfection:
The cuff can be completely inserted into liquid disinfectant. Only disinfectants with approved efficiency and in accordance with the national standards can be used. Before next use, please ensure that no liquid is remaining in the cuff. This can affect the measurement results negatively and it can damage the manometer technology.

IMPORTANT!
Do not iron nylon-velcro cuffs, Disinfectable one piece cuff. Never expose the cuffs to intensive solar radiation! Never touch the cuff covers or bladders with a sharp instrument, since this could cause damage!

Maintenance The product does not require any maintenance.

Test of precision
Remove the tube from the manometer and hold the manometer in a vertical position. When the pointer stands still at 0 on the scale, the instrument has been adjusted properly. If the pointer is below or above 0, the instrumt must be recalibrated. Either take it to an authorised **Riester** dealer or send it to us.

Monitoring of instruments
All countries except for Germany:

The respective legal provisions apply for all countries, except for Germany. The reference manometer, which is used for calibration, must be traceable to national and international measurement standards.

Warning: It is not allowed to make changes to the device.

Technical Data
Type of indication: Aneroid, desk, wall, stand and anaesthetic model
Round scale
Scale graduation: Increments of 2 mm Hg
Range shown on scale: 0 to 300 mm Hg
Measurement range: 0 to 300 mm Hg
Pointer movement per 1 mm deflation: 1 or 2, depending on the particular model
Tube adapter: Hand-, Tisch-, Wand-, Stand- und Anästhesiemodelle
Pressure generation: Bulb
Pressure reduction: Air-release valve that can be regulated

Storage conditions: 10°C (50°F) to 40°C (104°F) at a relative humidity of 85 % (non-condensing)
-20°C (4°F) to 70°C (158°F) at a relative humidity of 85 % (non-condensing)
Versions:
Type of indication: Aneroid, desk, wall, stand and anaesthetic model
Round scale
Scale graduation: Increments of 2 mm Hg
Range shown on scale: 0 to 300 mm Hg
Measurement range: 0 to 300 mm Hg
Pointer movement per 1 mm deflation: 1 or 2, depending on the particular model
Tube adapter: Hand-, Tisch-, Wand-, Stand- und Anästhesiemodelle
Pressure generation: Bulb
Pressure reduction: Air-release valve that can be regulated

Disposal: You can obtain information on disposal from your appropriate local facility or from your local environmental advisor.

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Avvertisements / Contre-indications:

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