SABAglyde®

INSTRUCTIONS FOR USE

EN English	DE Deutsch	CS Čeština
DA Dansk	EL Ελληνικά	ES Español
FI Suomi	FR Français	IT Italiano
LT Lietuvių	NL Nederlands	NO Norsk
PL Polski	PT Português	RO Română
RU Русский	SV Svenska	TR Türkçe
HU Magyar	UA українська	SL Slovenščina

SABANA®
MEDIZINBEDARF

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DESCRIPTION:
SABAglyde® suture (Poliglecaprone-25), is a sterile synthetic monofilament absorbable surgical suture prepared from a copolymer of glycolide and e-caprolactone. SABAglyde® sutures are available both in undyed and dyed variants. Undyed suture is beige in colour. Dyed variant of suture is violet in colour, which is dyed with D&C violet No 2 (conforming to U.S.Code of Federal regulations). SABAglyde® suture has been found to be non-antigenic, non-pyrogenic and elicits only a slight tissue reaction during the process of absorption. The composition of dyed variant of Poliglecaprone-25 suture is ≥ 99.93% of glycolide and e-caprolactone copolymer(75:25), and ≤ 0.07% of D&C violet No 2. The composition of undyed variant of Poliglecaprone-25 suture is 100% glycolide and e-caprolactone copolymer(75:25). SABAglyde® sutures comply with the requirements of the European Pharmacopoeia (E.P) for Sterile Synthetic Absorbable Monofilament Sutures, and United States Pharmacopeia (U.S.P) for Absorbable Surgical suture (Synthetic), except for slight over size in diameter.

INTENDED USE
SABAglyde® suture is intended for use in general soft tissue approximation and/or ligation where an absorbable suture material is indicated. The safety and effectiveness of SABAglyde® suture has not been established in cardiovascular tissue, neurological tissue, microsurgery and ophthalmic surgery.

INDICATION
Indicated in surgical procedures where a combination of an absorbable suture and short term wound support is required for wound closure, soft tissue approximation and/or ligation.

CONTRA INDICATIONS
SABAglyde® suture being absorbable, should not be used where extended approximation of tissue under stress is required. For example: fascia.

INTENDED USERS
Trained and registered health care professional only.

INTENDED PATIENT POPULATION
This suture can be used in patients irrespective of age and sex, in line with the intended use, indications and contraindications.

PERFORMANCE
SABAglyde® suture after implantation elicits a minimal initial inflammatory reaction in tissues and is eventually replaced with an ingrowth of fibrous connective tissue. SABAglyde® suture being absorbable in nature gets absorbed in the tissues by means of hydrolysis, during which the Poliglecaprone-25 polymer is degraded into monomers of adipic acid which is eventually absorbed and metabolized in the body. This process of absorption

initially leads to progressive loss of tensile strength, followed by a loss of mass.
Implantation studies in rats show the following profile:

Days post Implantation	**Approximate % Original Strength Remaining**
7 days	68% - 79%
14 days	39% - 41%

Absorption of SABAglyde® suture is essentially complete approximately in 90 days.

APPLICATION

Suture should be selected and implanted depending on patient's condition, surgical experience, surgical technique and wound size.

WARNINGS

- This suture is intended to be used by a trained and registered health care professional.
- Users should be familiar with surgical procedures and surgical suturing techniques involving absorbable sutures before using SABAglyde® suture for wound closure, as risk of wound dehiscence may vary with the site of application and the suture material used. Users should consider the in vivo performance when selecting a suture.
- Users must apply their professional judgment when determining the appropriate sized suture based on the specific indication, wound size, preferred surgical technique, patient's condition and history.
- The use of this suture may be inappropriate in elderly, malnourished, or debilitated patients, or in patients suffering from conditions which may delay wound healing.
- As with any foreign body, prolonged contact of this suture with salt solutions such as those found in the urinary or biliary tracts may result in calculus formation.
- As an absorbable suture, it may transiently act as a foreign body and as with any foreign body, presence of bacterial contamination may enhance bacterial infectivity. Thus, acceptable surgical practice should be followed with respect to drainage and closure of infected or contaminated wounds.
- Wound dehiscence may occur when suture fails to provide adequate wound support in closure of the sites where expansion, stretching or distension occur. Hence use of supplemental non-absorbable suture should be considered by the user in the closure of sites which may undergo expansion, stretching, or distension, or which may require additional support.
- Under certain circumstances, notably orthopaedic procedures, immobilization by external support may be employed at the discretion of the user.
- Consideration should be taken in the use of absorbable sutures in tissues with poor blood supply as suture extrusion and delayed absorption may occur.
- Conjunctival, cuticular and vaginal epithelial sutures remaining in place for extended periods may be associated with localized irritation and should be snipped off or removed as indicated.
- Subcuticular sutures should be placed as deeply as possible to minimize the erythema and induration normally associated with the absorption process.

PRECAUTIONS

- Care should be taken to avoid damage when handling surgical needles. Grasp the needle in an area one third (1/3) to one half (1/2) of the distance from the attachment end to the point. Grasping in the point area could impair the penetration performance & cause fracture of the needle. Grasping at the butt or attachment end could cause bending or breakage of the needle. Reshaping needles may cause them to loose strength and make less resistant to bending and breaking. Broken needles may result in extended or additional surgeries or residual foreign bodies.
- Users should exercise caution when handling surgical needles to avoid inadvertent needle stick injury. Discard used needles in 'Sharps' containers.
- When handling this suture, care should be taken to avoid damage from handling. Avoid crushing or crimping damage due to application of surgical instruments such as forceps or needle holders.

- When required, straighten sutures with a gentle pull, never crush or rub them.
- As with any suture material, adequate knot security requires the accepted surgical techniques of flat and square ties, with additional throws as warranted by surgical circumstances and the experience of the user. The use of additional throws may be particularly appropriate when knotting monofilament suture.
- Do not tie bands around pouches/packs of sutures as it may damage the pouches/pack.
- Do not re-sterilize, do not reuse. Reuse of this device (or portions of this device) may create a risk of product degradation, which may result in device failure and/or cross-contamination, which may lead to infection or transmission of blood-borne pathogens to patients and users.

SIDE EFFECTS

Side effects associated with the use of SABAglyde® suture include allergic response in certain patients, transient local irritation at the wound site, transient inflammatory foreign body response, mild inflammatory tissue reaction, erythema and induration during the absorption process of subcuticular sutures, and rarely wound dehiscence. Like all foreign bodies, SABAglyde® suture may enhance an existing infection.

CLINICAL BENEFITS:

The clinical benefits of this suture are, it being a monofilament suture allows easy tissue passage and ensures minimal/no tissue drag and trauma. It provides secure knot holding capacity, and being monofilament in nature doesn't have interstices in its structure leading to reduced incidence of bacterial colonization, when compared to braided absorbable sutures. When used for tissue approximation and/or ligation following injury or surgery, it helps in closing dead space, minimizes blood loss and infection, and also aid in wound healing. It provides adequate initial tensile strength as per the U.S.P and E.P requirements. It provides adequate short term wound support during the critical wound healing period, even with smaller diameter sutures. It has a predictable absorption rate enabling decision making regarding the application of sutures in relevant tissues and conditions, as against natural sutures like catgut which doesn't have a predictable absorption rate. It is absorbed by the process of hydrolysis enabling mild tissue reaction and inflammatory response. It leads to less reactive scars with a lower tendency towards hypertrophic scar formation, and hence provides good cosmetic results, and therefore very beneficial for closing skin incisions to provide a superior cosmetic outcome. This suture being absorbable in nature, eliminates the discomfort of removal of sutures, unlike non-absorbable sutures, thereby saving the surgeon's time and lessening patient anxiety and discomfort.

STERILITY

SABAglyde® suture is sterilized by ethylene oxide gas. Do not re-sterilize, do not reuse. Do not use if package is opened or damaged. Discard opened, unused sutures.

SUPPLY

SABAglyde® suture is available in various U.S.P. sizes and E.P Metric sizes. Sutures are supplied sterile in pre-cut lengths, both non-needled and attached to various needle type, shape and length. Sutures are packed in a printed box with quantity indicated on the box label.

STORAGE

Recommended storage conditions: Below 30° C, away from moisture and direct sunlight. Do not use after expiry date.

DISPOSAL

Discard used sutures and needles contaminated with blood and tissue in the container meant for infectious waste. Unused expired pouches should be incinerated or disposal should be done as per local regulations.

SYMBOLS USED ON LABEL:

Do not re-use	Batch code
Mfg. Date & Country of manufacture	Use-by date
Manufacturer	Medical Device
Poly(glycolide-co-caprolacton)	Sterilized using Ethylene oxide
Single sterile barrier system with protective packaging inside	Unique device identifier
GS 1 Data matrix barcode	Catalogue number

Upper limit of temperature	Caution
Keep away from sunlight	Do not resterilize
Keep dry	Consult instructions for use
Do not use if package is damaged	Authorized representative in the European Community
UK Responsible person	Importer
Distributor	

CE 2265 European conformity mark and Identification of Notified Body. The product meets the safety and performance requirements of EU Medical devices Regulation (EU) 2017/745

BASIC UDI-DI : 8903837H0101010103ZG

Summary of safety and clinical performance(SSCP) is available in the European database on medical devices (Eudamed), where it is linked to the BASIC UDI-DI : 8903837H0101010103ZG
https://ec.europa.eu/tools/eudamed

Adverse events should be reported to Healthium Medtech Ltd, or an in-country representative, and to the health authority where the event occurred. Email id: care@healthiummedtech.com

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Healthium Medtech Limited
Sabana Medizinbedarf GmbH
EC REP MED DEVICES LIFESCIENCES B.V.