SABANA **SILK**

SABANA®
MEDIZINBEDARF

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(EN) Instructions for Use

Non-absorbable surgical suture with or without needle

1. DESCRIPTION

SABANA Silk sutures are braided, non-absorbable, sterile surgical suture composed of an organic protein called fibroin. This protein is derived from the domesticated species Bombyx mori (B.mori) of the family Bombycidae. Silk for braided material is processed to remove the natural waxes and gums. Braided silk is coated with silicone and is available undyed (Natural Ivory) and dyed in black with Sulphol Black 1 (EP Vol III) (Color Index Number 53185) or Logwood Black (Color Index Number 75290). Logwood Black is a natural Tinctorial Wood Extract obtained from Haematoxylon Campechianum. Conforms to US Code of Federal Regulations 21 CFR 73.1410.
For virgin silk the sericin gum is not removed and holds the twisted filaments together. Virgin silk is available dyed blue with methylene blue (Color Index Number 52015) and dyed black with Logwood Black (Color Index Number 75290).
SABANA Silk suture is available in a range of gauge sizes and lengths attached to stainless steel needles of varying types and sizes.
SABANA Silk suture complies with the requirements of the European Pharmacopoeia for Sterile Non-absorbable Silk Sutures and the United States Pharmacopoeia monograph for Non-absorbable Sutures.

2. INDICATIONS

SABANA Silk suture is indicated for use in general soft tissue approximation and/or ligation, including use in ophthalmic procedures.

Application

SABANA Silk sutures should be selected and implanted depending on patient condition, surgical experience, surgical technique and wound size.

Material performance

SABANA Silk suture elicits an initial inflammatory reaction in tissues, which is followed by gradual encapsulation of the suture by fibrous connective tissues. While silk is not absorbed, progressive degradation of the proteinaceous silk fibre in vivo may result in gradual loss of all the suture's tensile strength over time.

3. CONTRAINDICATIONS

The use of this suture in contraindicated in patients with known sensitivities or allergies to silk. Due to the gradual loss of tensile strength which may occur over prolonged periods in vivo, silk suture should not be used where permanent retention of tensile strength is required. It`s not indicated for the central circulatory system and central nervous system.

4. WARNINGS/PRECAUTIONS/INTERACTIONS

- Do not use if package is opened or damaged. Discard opened unused sutures.
- Do not use after expiry date.
- Users should be familiar with surgical procedures and techniques involving non-absorbable sutures before employing SABANA Silk sutures for wound closure, as risk of wound dehiscence, may vary with the site of application and the suture material used.
- Acceptable surgical practice should be followed for the management of infected or contaminated wounds.
- As with any exogenous material, prolonged contact of any suture with salt solutions, such as those found in the urinary or biliary tracts may result in calculus formation.
- Adequate knot security requires the standard surgical technique of flat and square ties with additional throws as indicated by surgical circumstance and the experience of the surgeon.
- When handling this or any other suture, care should be taken to avoid damage. Avoid crushing or crimping damage due to the application of surgical instruments such as forceps or needle holders.
- Care should be taken to avoid damage when handling surgical needles. Grasp the needle in an area one-third (⅓) to one half (½) of the distance from the attachment end to the point. Grasping in the point area could impair the penetration performance and cause fracture of the needle. Grasping at the butt or attachment end could cause bending or breakage. Reshaping needles may cause them to lose strength and be less resistant to bending and breaking.
- Users should exercise caution when handling surgical needles to avoid inadvertent needle stick injury.
- Discard used needles in "Sharps" container.
- Dispose of material in accordance with all the state, local, and hospital regulations. Responsibility for proper waste disposal is with the owner of the waste.
- Do not reuse: Infection hazard for patients and/or users and impairment of products functionality due to reuse. Risk of injury, illness or death due to contamination and/ or impaired functionality of the product.
- Do not resterilize: Infection hazard for patients and/or users and impairment of products functionality due to use of resterilized suture. Risk of injury, illness or death due to contamination and/or impaired functionality of the product.

5. ADVERSE REACTIONS

Adverse effects associated with the use of SABANA Silk include allergic response in patient that are known to be sensitive to silk, initial inflammatory tissue reaction, and transient local irritation at the wound site. Like all foreign bodies silk may enhance an existing infection.

6. STERILITY

SABANA Silk sutures are sterilized by ethylene oxide gas.
Do not resterilize!
Do not use if package is opened or damaged! Discard opened unused sutures!

7. STORAGE CONDITIONS

Recommended storage conditions for SABANA Silk: store at a temperature below 25°C, away from moisture, corrosion and direct heat.
Do not use after expiry date.

8. LABELING SYMBOLS

- yyyy-mm Date of manufacture (year-month)
- Dyed, woven and non-absorbable suture material
- Dyed, woven, coated and non-absorbable suture material
- Undyed, woven, coated and non-absorbable suture material
- LOT Batch designation
- STERILE EO Sterilized with ethylene oxide
- Do not reuse
- Do not resterilize
- Comply with the instructions for use
- Caution, consult accompanying document
- Keep dry
- yyyy-mm Use by date (year-month)
- Store in a dark location away from sunlight
- 25°C Upper temperature limit for storage
- Do not use if the packaging is damaged
- CE 0123 CE mark with code number of the Notified Body (TÜV SÜD Product Service GmbH)

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CE 0123 Oznaka CE z identifikacijsko številko priglašenega organa (TÜV SÜD Product Service GmbH)