SABApol®

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EN Instructions for Use

Sterile, synthetic, non-absorbable surgical suture material

1. DESCRIPTION

SABApol® is a monofilament, synthetic, non-absorbable, sterile surgical suture composed of an isotactic crystalline stereoisomer of polypropylene, a synthetic linear polyolefin. The molecular formula is $(C_3H_6)n$. SABApol® suture is available undyed (clear) and dyed blue with phthalocyanine blue(Color Index Number 74160).
SABApol® suture is available in a range of gauge sizes and lengths attached to stainless steel needles of various types and sizes.
SABApol® suture complies with the requirements of the European Pharmacopoeia for Sterile Non Absorbable Polypropylene suture and the requirements of the United States Pharmacopoeia monograph for Non absorbable Sutures.

2. INDICATIONS

SABApol® suture is indicated for use in general soft tissue approximation and/or ligation, including use in cardiovascular, ophthalmic and neurological procedures.

Application
Sutures should be selected and implanted depending on patient condition, surgical experience, surgical technique and wound size.

Material Performance
SABApol® suture elicits a minimal initial inflammatory reaction in tissues which is followed by gradual encapsulation of the suture by fibrous connective tissue. It is a non-absorbed suture. Due to its relative biological inertness it is recommended for use where the least possible tissue reaction is desired. As a monofilament it has been successfully employed in surgical wounds which subsequently become infected or contaminated where it can minimize later sinus formation and suture extrusion. Because of its lack of adherence to tissue, Polypropylene suture is effective as a pull-out suture.

3. CONTRAINDICATIONS

None known.

4. WARNINGS / PRECAUTIONS / INTERACTIONS

- Users should be familiar with surgical procedures and techniques involving non-absorbable sutures before employing polypropylene suture for wound closure, as risk of wound dehiscence, may vary with the site of application and the suture material used.
- As with any foreign body, prolonged contact of any suture with salt solutions, such as those found in the urinary or biliary tracts may result in calculus formation.
- Acceptable surgical practice should be followed for the management of infected or contaminated wounds.
- Care should be taken to avoid damaging the surface of the material with surgical instruments as this could lead to fracture of the material in use. Avoid crushing or crimping damage due to the application of surgical instruments such as forceps of needle holders.
- Adequate knot security requires the standard surgical technique of flat and square ties with additional throws as indicated by surgical circumstance and the experience of the surgeon.
- The use of additional throws is particularly appropriate when knotting polypropylene sutures.
- Care should be taken to avoid damage when handling surgical needles. Grasp the needle in an area one-third (⅓) to one half (½) of the distance from the attachment end to the point. Grasping in the point area could impair the penetration performance and cause fracture of the needle. Grasping at the butt or attachment end could cause bending or breakage. Reshaping needles may cause them to lose strength and be less resistant to bending and breaking.
- Users should exercise caution when handling surgical needles to avoid inadvertent needle stick injury.
- Discard used needles in "Sharps" container.
- Dispose of material in accordance with all the state, local, and hospital regulations. Responsibility for proper waste disposal is with the owner of the waste.
- Do not reuse: Infection hazard for patients and/or users and impairment of products functionality due to reuse. Risk of injury, illness or death due to contamination and/or impaired functionality of the product.
- Do not resterilize: Infection hazard for patients and/or users and impairment of products functionality due to use of re-sterilized suture. Risk of injury, illness or death due to contamination and/or impaired functionality of the product.
- Do not use after expiry date.

5. ADVERSE REACTIONS

Adverse reactions associated with the use of SABApol® include minimal initial inflammatory tissue reaction and transient local irritation at the wound site. Like all foreign bodies Polypropylene suture may potentiate an existing infection.

6. STERILITY

SABApol® sutures are sterilized by ethylene oxide gas.
Do not re-sterilize!
Do not use if package is opened or damaged! Discard opened unused sutures!

7. STORAGE CONDITIONS

Recommended storage conditions for SABApol®: store at a temperature below 25°C, away from moisture, corrosion and direct heat.
Do not use after expiry date.

8. LABELING SYMBOLS

yyyy-mm Date of Manufacture (year, month)

Dyed, monofilament, non-absorbable suture material

Undyed, monofilament, non-absorbable suture material

LOT Batch designation

STERILE EO Sterilized with ethylene oxide

Comply with the instructions for use.

Caution, consult accompanying document

Do not reuse.

Do not resterilize.

yyyy-mm Use by date (year-month)

Store in a dark location away from sunlight.

Keep dry

25°C Upper temperature limit for storage

Do not use if the packaging is damaged.

CE 0123 CE mark with code number of the Notified Body

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