

# Pen options

NovoFine® needles fit all Novo Nordisk, and most diabetes pens.



FlexTouch®



NovoPen® 5



FlexPen®



NovoPen® 4



NovoPen Echo®

NovoFine® needles also fit most diabetes pens, such as:

- |                  |                    |                       |
|------------------|--------------------|-----------------------|
| KwikPen®         | • Autopen® 24      | • BerliPen® 302       |
| SoloStar®        | • Omnicar® Pen 31  | • HumaPen® Luxura™    |
| Liprolog® Pen    | • OptiPen® Pro 1   | • HumaPen® Luxura™ HD |
| Huminsulin® Pen  | • Ypsopen®         | • HumaPen® MEMOIR™    |
| Exenatide Pen    | • BerliPen® Areo   | • OptiClik®           |
| Autopen® Classic | • BerliPen® Areo 2 | • TactiPen®           |

Brands are trademarks of their respective owners.

## The NovoFine® range

Available in packs of 100

Length	Diameter
4 mm	32G*
6 mm	32G*
6 mm	31G*
8 mm	30G*



8 mm

30G\* Autocover®

\* The G stands for gauge and is a measure of needle thickness; the higher the number, the thinner the needle.