



TOP-NOTCH  
PERFORMER

**iFuge M15K**  
(Table Top Centrifuge)

# iFuge M15K

(Table Top Centrifuge)

In any research facility or lab where accuracy reigns supreme, you need an equipment that can deliver precision to the last decimal point. Such an equipment might be hard to get hold of but not impossible, especially when Neutation is here to help you.

With Neutation's iFuge M15K you can expect a lot as it delivers even more. The lightning fast 15000 RPM microprocessor controlled centrifuge can perform at peak efficiency delivering top-notch accuracy in every cycle despite high load.



## Witness The Champion



Brushless Dc Motor For  
Maintenance Free Long Life



Variable Speed  
From 500-15000 RPM



Microprocessor  
Controlled Efficient  
Performance.



Intuitive Interface For  
Quick And Convenient  
Setting With Last Run  
Memory Feature



One Touch RPM / RCF  
Conversion Operation

Variable Timer Setting  
From 1 To 99 Min &  
Infinite Mode

### Witness The Champion



Lid-lock Safety Feature  
Opens The Lid Automatically  
After Every Run



Enhanced Airflow Design  
For Less Noise And  
Greater Sample Safety



Smooth & Silent Operation  
Even At Maximum Speed

Safe Low  
Voltage Device



Compact Design And  
Can Also Be Used In  
Hoods And Cold Rooms



### ROTOR OPTION



12 x 1.5/2 ml | 15596 x g



2 x 8 x 0.2 ml strip rotor

# Specifications

| iFuge M15K   |   |
|--|---|
| Specifications                                     | Details   |
| Rotor Capacity                                     | 12 x 1.5ml / 2.0ml Microtubes   |
| Maximum RPM / RCF                                  | 15000 / 15596g  |
| Display  | Large Digital   |
| Imbalance Detection                                | Yes   |
| Motor Type   | (BLDC) Brushless DC motor   |
| Timer Setting                                      | Countdown Timer From 1-99 Mins & Infinity Mode  |
| Dimensions (WxDxH)                                 | 187 x 265 x 115 mm  |
| Noise Level  | < 60 db   |
| Weight   | 2.2 Kg  |
| Rotors and Adaptors<br>(included in standard pack) | Reduction Adaptors For 0.2 ml Microtubes,<br>Reduction Adaptors For 0.5 ml Microtubes |

# Ordering Information

| iFuge M15K   |                       |   |
|--------------|-----------------------|---|
| Product Name | Cat #                 | Details                                       |
| iFuge M15K   | M112-5115 UK/AU/US/EU | High Speed Micro Centrifuge<br>Upto 15000 RPM |

| Accessories             |              |  |
|-------------------------|--------------|--|
| Part Name               | Cat #        | Details                                |
| 12 Slot Rotor           | M12-R12      | 12 x 1.5/ 2 ml tubes slot rotor        |
| Rotor Lid               | M12-R12E     | Rotor Lid for 12 slot rotor            |
| Rotor Nut               | M12-RN       | Rotor Nut for the Rotor Lid            |
| PCR Strip Rotor         | M12 - R2 x 8 | 2 PCR stirp rotor in Black colour      |
| PCR Strip Rotor Lid     | M12-R2 x 8E  | -                                      |
| Adaptors for 0.2 ml     | M12-D1       | Reduction Adaptors of 0.2 ml tubes     |
| Adaptors for 0.4/0.5 ml | M12-D2       | Reduction Adaptors of 0.4/0.5 ml tubes |
| Power Supply adaptor    | M12-A12      | Power adaptor with specified pin type  |

**NEUATION**  
TECHNOLOGIES  
www.neuation.com

AN ISO 9001 : 2015 & ISO 13485 : 2016 CERTIFIED MFG. COMPANY

\*We also execute OEM projects. Specifications can be changed without notice for quality improvement.