

DERMLITE DL5 TECHNICAL DATA AND COMPARISON TABLE

	DL1	DL100	DL200 HR	DL200 Hybrid	DL3N	DL4W	DL4	DL5
POLARIZATION	Polarized & Non-Pol.	Polarized	Polarized	Polarized & Non-Pol.	Polarized & Non-Pol.	Polarized & Non-Pol.	Polarized & Non-Pol.	Variable polarization (incl. polarized/non-polarized toggle)
OPTICS	Ø15mm 10x	Ø15mm 10x	Ø25mm 10x	Ø25mm 10x	Ø25mm 10x	Ø30mm 10x	Ø30mm 10x	Ø32mm 10x
ILLUMINATION	4 white LEDs (pol./non-pol.)	8 white pol. LEDs	21 white pol. LEDs	21 white LEDs (15 pol., 6 non-pol.)	18 white LEDs (12 pol., 6 non-pol.)	24 white LEDs (18 pol., 6 non-pol.)	18 white LEDs (12 pol., 6 non-pol.), 6 orange pol. LEDs	16x white, 4x UV (365nm), 4x Pigment Boost, 1 torch LED
INFECTION CONTROL	-	-	IceCap®	IceCap	-	IceCap	IceCap	IceCap
POWER	Rechargeable	Disposable (Type 2CR5)	Rechargeable	Rechargeable	Rechargeable, includes charging base	Rechargeable	Rechargeable	Rechargeable, includes base
SMARTPHONE USE	Yes	-	Yes	Yes	Yes	Yes	Yes	includes MCC adapter
CONTACT PLATE	Bayonet	-	Snap-on	Snap-on	Fixed	Snap-on	Snap-on	Bayonet style

INCLUDED ACCESSORIES	MCC smartphone adapter	-	Silicone Sleeve with Lanyard, IceCap supply, USB cable	Silicone Sleeve with Lanyard, IceCap supply, USB cable	Desktop Charger	USB Cable	USB Cable	Desktop charging base, phone adapter, magnetic eyepiece , 100 mm ruler, leather belt pouch, protective holster, 5 IceCaps, 2m USB-A to USB-C charging cable, cleaning cloth
					Leather Pouch	Belt Holster	Belt Holster	
					Eyepiece	Silicone Sleeve IceCap sample supply	Silicone microfiber Sleeve IceCap sample supply	