



Item number 1999918623



Classic

Visitor jacket, ABENA Classic, S/M, 100x72cm, white, PP, non-sterile, single use

- ✓ Soft and comfortable
- ✓ Long sleeves with elastic and hook-and-loop fastening
- ✓ Unisex design
- ✓ Minimizes cross-contamination
- ✓ Hygienic
- ✓ Space-saving



Product Description

Long sleeve visitor jacket with elastic sleeves and hook-and-loop fastening. Hook and loop tape is easy to pull apart and re-seal again. Used to protect the user and environments in which high hygienic standards are required.



Specifications

Base name	Visitor jacket
Brand	ABENA
Sub-brand	Classic
Size	S/M
Color	White
Features	Non-sterile
Single or multiple use	Single use
Material	PP
Length/depth	100 cm
Width	72 cm
Thickness	25 g/m2
Certifications	CE. Food contact materials. CAT I.
CE Category (Personal Protective Equipment)	CAT I
Directives, regulations and acts	(EU) No 10/2011, (EC) No 1935/2004, (EC) No 2023/2006, (EU) 2016/425, BEK nr 681 af 25/05/2020
Because of the products composition it has an expiry date from the production date on:	5 years
Storage Instructions	Store dry, room temperature and no direct sunlight.
Product Disposal Instructions	Can be disposed of with normal household waste sorted according to local regulations.
Packaging Disposal Instructions	Can be recycled or incinerated.

Packaging data

Unit	Contains	Length	Width	Height	EAN
cll	10 pck	45.5 cm	33.6 cm	34.6 cm	5703538496516
pck	10 pcs	32.8 cm	32.4 cm	4.2 cm	5703538496509



The CE mark guarantees that a product is safe to use and complies with all safety precautions. CE stands for Conformité Européenne (European Conformity) and is mainly found on electronic equipment, safety equipment, construction products, and medical devices.



CAT I PPE is only for products with minimal risk.



The glass-fork-symbol guarantees that products live up to the demands stated in the European legislation for food contact materials. The symbol is not mandatory but advised on products used for food contact.



The product does not contain latex.



The product is liquid repellent.



The product is only for single-use.