

Askina® Skin Freeze

CRYOTHERAPEUTIC KIT FOR PROFESSIONAL USE | TECHNICAL DATA



PRODUCT CLASSIFICATION

Active medical device, CE class IIa

QUALITY CERTIFICATES

EN ISO 13485:2016 (CE 0344)

MANUFACTURER

 Koninklijke Utermöhlen NV
De Overweg 1, 8471 ZA Wolvega, The Netherlands

DISTRIBUTOR

 B. Braun Melsungen AG
Carl-Braun-Str. 1, 34212 Melsungen
Germany

Askina® Skin Freeze

CRYOTHERAPEUTIC KIT FOR PROFESSIONAL USE | TECHNICAL DATA

PORTABLE CRYOSURGICAL SYSTEM FOR THE REMOVAL OF WARTS AND OTHER SKIN LESIONS

KIT COMPONENTS	aerosol can filled with 170 ml cryogenic liquid box with applicators instructions for use REF 9380702 9380703 9380704 packaging shelf		
DEVICE VARIANTS	life storage conditions	Applicator sizes	Applicators/pack
		2 mm tip for small lesions	60
		5 mm tip for standard therapy	50
		Mix of 2 and 5 mm applicators	2 x 30
LOGISTICS	kit package for 50–60 treatments, carton with 6 kits 3 years		
			
MATERIAL DATA		170 ml aerosol can	
		can	aluminum
		spraying device	polyethylene
	 For hazard and safety information: please request safety data sheet.	Cryogenic liquid DANGER	dimethyl-ether, propane, butane and iso-butane Extremely flammable aerosol. Pressurized container: May burst if heated. Keep out of reach of children. Do not breathe spray. Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. Do not pierce or burn, even after use. Do not spray on an open flame or other ignition source. Protect from sunlight. Do not expose to temperatures exceeding 50°C. Dispose of contents/container in accordance with local regulations.
			free of natural rubber latex
		applicator stick	polypropylene
		tip	cotton bead, polyurethane foam cover