


# Typhoid Antigen Rapid Test Kit (Colloidal Gold Method)

The Typhoid Antigen Rapid Test Kit is a lateral flow immunoassay assay designed for qualitative detection of Salmonella Typhi (S. typhi) antigens in human fecal, serum or plasma. This method is easy to operate, and the result can be obtained by observing whether there is a red and purple line in 15-20 minutes. It is intended to be used as a screening test and provides a preliminary test result to aid in the diagnosis of infection with S. typhi.


## Test Procedure

- Specimen: Fecal, serum or plasma.
- Specimen collection: Specimens of all types are collected by conventional methods.
- Test method:

### Test method:



### Interpretation of test result



## Specifications

Components	Test card/ Sample extraction tubes /Sample diluent/Sample droppers
Pack Size	25Tests/Kit; 50 Tests/Kit
Storage Conditions	2°C ~ 30°C ;Valid for 24 months
Test samples	Fecal, Serum & Plasma
Reaction conditions	15~30°C 15~20 minutes
Special instrument	No need

## Product Informations

Product Name	Cat.No	Package	Sample Diluent	Sample Droppers	Sample Type
Typhoid Antigen Rapid Test Kit (Colloidal Gold Method)	BSK15S1S	25 Tests/Kit	1 × 3 mL/Tube	25	Fecal, Serum & Plasma
	BSK15M1S	50 Tests/Kit	2 × 3 mL/Tube	30	