


# Typhoid IgG/IgM Rapid Test Kit (Colloidal Gold Method)

The Typhoid IgG/IgM Rapid Test Kit is a lateral flow immunoassay assay designed for qualitative detection and differentiation of IgG and IgM Typhoid antibodies produced against Salmonella Typhi (S. typhi) in human whole blood, serum or plasma. This easy-to-use format with rapid results allows point-of-care testing, providing information to assist with treatment and patient management decisions.


## Test Procedure

- Specimen: Whole blood, serum or plasma.
- Specimen collection: Specimens of all types are collected by conventional methods.
- Test method:

Test method:



Interpretation of test result



## Specifications

Components	Test card/Sample diluent/ Sample Droppers
Pack Size	25Tests/Kit; 50 Tests/Kit
Storage Conditions	2°C ~ 30°C ;Valid for 24 months
Test samples	Whole blood, Serum or Plasma
Reaction conditions	15~30°C 15~20 minutes
Special instrument	No need

## Product Informations

Product Name	Cat.No	Package	Sample Diluent	Sample Droppers	Sample Type
Typhoid IgG/IgM Rapid Test Kit (Colloidal Gold Method)	BSK11S1S	25 Tests/Kit	1 × 3 mL/Tube	25	Whole blood, serum & plasma
	BSK11M1S	50 Tests/Kit	2 × 3 mL/Tube	50	