



germanengineering

PO 80 Pulsoximeter

Ermittlung der arteriellen Sauerstoffsättigung (SpO) und der

Herzfrequenz (Pulsfrequenz)	2
Sehr einfache und vollkommen schmerzfreie Messung	
Klein und leicht für zu Hause und unterwegs	
Besonders geeignet für Personen mit:	
• Herzinsuffizienz • Chronisch obstruktiven Lungenerkrankungen • Asthma bronchiale • Schlafapnoe	
Anwendbar bei Sport in großen Höhen (z.B. Bergsteigen, Skifahren und Sportfliegerei)	
Alarmfunktion - Grenzwerte individuell einstellbar	
Kontinuierliche Aufzeichnung bis zu 24h	
PC-Software zur detaillierten Darstellung und Auswertung	
Betriebssysteme: - Win XP - Win Vista - Win 7	
USB-Kabel und USB-Ladegerät	
Aufladbar per USB (integrierter Lithium-Ionen-Akku)	
XL Farbdisplay mit 4 Ansichtsperspektiven	
Einstellbare Display-Helligkeit	
Graphische Pulsanzeige und zuschaltbarer Pulston	
Batteriestandsanzeige	
Abschaltautomatik	
Inkl. Halteband und Gürtel-Tasche	
Medizinprodukt	
Produktmaße: 57 x 32 x 30 mm	
Produktgewicht: ca. 42 g	
3 Jahre Garantie	
VE: 4 / Transportkarton: 8 x 4	
EAN-Nr.: 421112545440 1	
Art.-Nr.: 454.40	

PO 80 Pulse Oximeter

Determining arterial oxygen saturation (SpO) and

heart rate (pulse frequency)	2
Simple and completely pain-free measurement	
Small and light – for use at home and on the move	
Particularly suitable for persons with:	
• Heart failure • Chronic obstructive pulmonary diseases • Bronchial asthma • Sleep apnoea	
Suitable for sports at high altitudes (e.g. mountaineering, skiing and aviation sports)	
Alarm function - limit values can be individually set	
Continuous recording of up to 24h	
PC software for detailed display and evaluation	
Operating system: - Win XP - Win Vista -Win 7	
USB cable and USB charger	
Can be charged by USB (integrated lithium-ion battery)	
XL colour display with 4 available views	
Adjustable display brightness	
Graphical display of pulse and optional pulse tone	
Battery status indicator	
Automatic switch-off	
Incl. retaining strap and belt pouch	
Medical product	
Product measurements: 57 x 32 x 30 mm	
Product weight: approx. 42 g	
3 year guarantee	
Sales unit: 4 / Shipping carton: 8 x 4	
EAN no.: 421112545440 1	
Item no.: 454.40	