

Specifications : BT-410F Fixed Beamlight Headlamp

Functional Characteristics

LED Light

Illumination (at 25cm)	15,000 ~ 30,000lx
Color temperature	6,000Kelvin
Life time	50,000hrs

Function

Battery state display	Green / red(low battery)
Battery module attachable	
Light Intensity Adjustable	

Power

Adaptor	Input : AC 100 ~ 240V (50/60Hz)
	Output : DC 4.2V/750mA 3.7V,
Rechargeable bBattery	2,600mA Li-ion Operating time :
	4hrs Charging time : 4hrs

Standard Configurations

Operation manual	1ea
Battery module	1ea
Adapter	1ea

Optiopns

Loupe	x 3.5 (diopter)
Battery module	

Physical Characteristics

Dimension

Headband	534 ~ 638mm 282(W) x
Packing	187(D) x 117(H)mm 600(W) x
Carton box (5ea)	300(D) x 215(H)mm

Weight

Headband	119g
Battery module	69g
Packing	670g
Carton box (5ea)	4.35Kg

Environmental Conditions

Operating temperature	10 ~ 40°C (50 ~ 104°F)
Operating humidity	5 ~ 85% non-condensing
Storage temperature	-20 ~ 60°C (-4 ~ 140°F)
Storage humidity	0 ~ 95% non-condensing

Warranty

1year	
-------	--

Certificates

CE, FDA, CE, FDA	
------------------	--