


CETRO® Anoscope

1542




CETRO® Anoscope

Anoskop 1542 for single use is made in plastic. The anoscope is used in colorectal examinations and is constructed of a material that conducts light very well and provides an optimal view for the examiner. The anoscope's carefully rounded edges and surfaces contribute to a gentle examination for the patient.

- The beveled end provides better access to hemorrhoid treatments.
- Fits with the same light fixture as for rectoscopes and tube proctoscopes.
- Fits with several of the lighting systems on the market today.
- Melting temperature: 100 °C
- Not autoclave able
- Latex free

Material: Polyethylene



ARTICLE NO.	PRODUCT	DEPT. PACKAGING	TRANSP. PACKAGING	TRANSP. PACKAGING
1542	8 g	25 pcs	150 pcs	1,6 kg

