

CETRO® Rectoscope

1561




CETRO® Rectoscope

Rectoscope 1561 for single use. The tube is constructed in a material that conducts light very well and provides an optimal view for the examiner. The rectoscope's smoothly rounded edges and surfaces contribute to a gentle examination for the patient. Compatible with Heine and Welch Allyn.

- Stable and easy-to-grip construction
 - Clear scale
 - Suitable for different light heads
 - Very small transition between obturator and tube
 - The light transmission provides correct tissue color and a clear view
- Melting temperature: 100 °C
 - Not autoclave able
 - Latex free

Materials Glass-clear polystyrene, Polyethylene HD, Polyethylene LD



ARTICLE NO.	PRODUCT	DEPT. PACKAGING	TRANSP. PACKAGING	TRANSP. PACKAGING
1561	50,5 g	25 pcs	150 pcs	10 kg

