


CETRO® Tube Proctoscope

1551




CETRO® Tube Proctoscope

Tube proctoscope for single use. The tube is constructed of a material that conducts light very well and provides an optimal visibility for the examiner. The tubular proctoscope's rounded edges and surfaces contribute to one gentle examination for the patient.

- Stable and grip-friendly construction
 - Clear scale
 - Fits different light heads
 - Very small transition between obturator and tube
 - The light transmission gives the correct tissue color and free visibility
- Melting temperature: 100 ° C
 - Not autoclavable
 - Latex Free

Material: Polyeten HD



ARTICLE NO.	PRODUCT	DEPT. PACKAGING	TRANSP. PACKAGING	TRANSP. PACKAGING
1551	32 g	25 pcs	150 pcs	5 kg

