

CETRO® Spatula

1375




CETRO® Spatula

Disposable spatula for Gynecological Sampling. A spatula that has an ideal length of 20 cm for best handling, stable and comfortable to hold.

- Good adapted length
- Melting Temperature: 100 °C
 - Not autoclaveable
 - Latex free

Material: Polypropylene



ARTICLE NO.	PRODUCT	DEPT.PACKAGING	TRANSP.PACKAGING	TRANSP.PACKAGING
1375	3 g	100 pcs	1200 pcs	5 kg

