



COMeN

COMeN Share with the World



STAR5000E

Fetal & Maternal Monitor



Shenzhen Comen Medical Instruments Co., Ltd.
Add: No.2 of FIYTA Timepiece Building, Nanhuan Avenue, Matian Sub-district,
Guangming District, Shenzhen, 518106, Guangdong, China.
Tel:+86-755-2643 1236 400-700-9488 Fax:+86-755-2643 1232
Website: en.comen.com E-mail: info@szcomen.com
P/N: EN-STAR5000E-4P-20200302-V1.3

STAR5000E Fetal & Maternal Monitor



- Folding-up screen, easy for observation.



- 5.6" LED screen.
- Built-in high-capacity Li-ion battery, above 2 hours working time.
- 150mm wide printer. Special extractable paper slot, easy for paper replacement and observation.
- Support 3cm/Min, 2cm/Min, 1cm/Min
- real-time printing.
- Support timing printing.
- 3-level audio/visual alarms.
-



Support data storage and software upgrading by external USB disk. External SD card optional, for data transferring.



Portable hidden handle.

- Advanced patient files management system.
- Manual/automatic FM counting, automatic fetal status scoring system.
- FHR signal strength indicator, twins cross channel verification.
- 60 hours real-time waveform review
- Timing monitoring function, avoid overtime monitoring.
- Wire & wireless central monitoring station connection.
-

