



# High-Definition Wireless Ultrasound for **OB/GYN**

## Streamline Workflows

With best-in-class imaging and affordable pricing, Clarius handheld ultrasound is the ideal system to add to your OB/GYN practice for quick exams to reassure parents to be.



# Improve Diagnostic Confidence with AI-Powered Imaging

Why wait for access to a bulky ultrasound system? Affordable and wireless, Clarius streamlines OB/GYN exams. Ideal for first visits, quick check-ups, and when you need to investigate pelvic pain, pre-term labor symptoms, and fetal viability, the Clarius C3 HD3 convex and EC7 HD3 endocavity scanners deliver the image quality you need for a confident assessment.

- Confirm fetal viability and location of the pregnancy
- Detect early intrauterine pregnancy under 10 weeks
- Estimate gestational age and measure heart rate
- Assess fetal lie, placental position in third trimester
- Perform biometric measurements for singleton or multiple pregnancies
- Determine amniotic fluid volume and Biophysical Profile scoring
- Confirm IUD position
- Perform antral follicle count for IVF




**Clarius C3 HD3  
Curvilinear Scanner**

**Frequency:** 2-6 MHz  
**Max Depth:** 40 cm




**Clarius EC7 HD3  
Endocavity Scanner**


**Frequency:** 3-10 MHz  
**Max Depth:** 15 cm




Fetal Spine



Ovary (C3)




Ovary (EC7)



Fetal AC



SAX Kidney



Diaphragm



I find that I'm using it more often because it seamlessly fits in as an extension of the physical exam. I'm very impressed. It's quite a bit smaller and lighter, so it fits very nicely into my hand and is very simple and easy to use. It's very sleek. The image quality's been excellent. It's very crisp and clean and definitely in comparison to the previous cart-based ultrasound we were using, it's been night and day."

**Dr. Ulrike Dehaeck, MD**

# High-Resolution Ultrasound in a Portable Wireless Scanner



## High-Definition

Octal Beamforming technology and 192 piezoelectric elements deliver high performance and resolution.



## Easy To Use

Friendly AI-powered app for your Apple or Android automatically delivers the best imaging.



## Wireless

Steadfast connection. Total freedom. Without cords, cleaning and disinfection are simple and fast.



## Specialized

AI-powered with specialized workflows for obstetrics.



## Ultra-Portable

Ideal for small spaces, Clarius fits in your pocket. Take it anywhere!




## Clarius Cloud

Manage your exams and build reports anywhere, anytime with unlimited Clarius Cloud storage.



The image quality of the Clarius on our iPad is far superior to what we had before. It is easy to enlarge images to point things out to patients. The touch screen is fantastic and very intuitive. This is a great product."

**Dr. Brian Driedger, MD**



# Advanced OB Package

Get the additional flexibility you need with additional tools for more advanced obstetrical scanning. The new labor and delivery preset optimizes imaging during labor with tools like angle of progression.




## Advanced Measurements

Measure and report on gestational age through packages for long bones, cranial, and trunk anatomy.



## Growth Tables

Automatically generate comprehensive growth charts included within PDF reports.



## Act One™

Delight patients by emailing them an animated movie automatically created from select images and video clips.

## Additional Features

- Advanced Outputs
- DICOM
- Elastography
- Pulsed Wave Doppler
- Procedure Recording Mode
- Specialized Presets
- Split Screen Mode

# Clarius Membership

The Advanced OB Package is included with your Clarius Membership.

## Membership Bundle

With your membership plan, you'll receive a discount on a scanner and the following:

- ✓ Support
- ✓ 60 Minutes of 1:1 Training
- ✓ Unlimited Clarius Cloud Storage
- ✓ Clarius Live
- ✓ Clarius Education
- ✓ Advanced OB Package

## One-Time Purchase

You may also purchase your scanner outright with the Advanced OB Package à la carte:

- ✓ Support
- ✓ 60 Minutes of 1:1 Training
- \$ Advanced OB Package

[REQUEST A QUOTE](#)

Clarius is on a mission to make accurate, easy-to-use, affordable ultrasound available to all medical professionals in every specialty. Our ultra-portable wireless scanners deliver best-in-class imaging with an easy-to-use app for iOS and Android that is powered by artificial intelligence and connected to the cloud. More than two million high-definition scans have been performed using Clarius wireless scanners, now available in over 90 countries worldwide.

## Clarius Mobile Health

2985 Virtual Way #130,  
Vancouver, BC, V5M 4X7  
Canada  
hello@clarius.com  
Phone 1-778-800-9975

MKTG-00217 Rev 1