



High-Definition Wireless Ultrasound for **Facial Aesthetics**

Improve Patient Safety

Clarius handheld ultrasound is the leading choice for aesthetics practitioners to clearly visualize facial and superficial anatomy in real-time to safely guide procedures, including cosmetic fillers.




Image Clarity You Can Trust for Facial Aesthetics

Medical aesthetic practitioners looking to improve patient safety rely on the Clarius L20 HD –
3

the world's only ultra-high frequency wireless ultrasound with an Advanced Aesthetic Package. Facial vascular mapping takes just minutes. Impress your patients with clear imaging of the skin, vessels, muscles, and fascia for safe, consistent filler treatments. You'll gain confidence by visualizing your needle placement to avoid vascular occlusions and dissolve unwanted fillers.

- Perform facial mapping to avoid vascular occlusions
- Confirm filler placement and vascularity after injections
- Visualize your needle and nodules to treat complications like vascular occlusions using ultrasound guidance
- Measure dermal thickness prior to microneedling
- Visualize vessels for reconstructive/ cosmetic facial surgery




Clarius L20 HD3
Ultra-High Frequency

Frequency: 8-20 MHz


Max Depth: 4 cm

Weight: 290g (0.64lb)


Dimensions: 147 x 76 x 32mm (5.8 x 3.0 x 1.3in)




Labiomental Artery




Facial Artery




Filler Around Eye



Temporall Artery



Lab Artery



Temporal Vessels



When I started using the Clarius L20 HD, I saw firsthand the benefits of using ultrasound every day. It's better for my patients because I can place fillers exactly where I want and it's better for my peace of mind because I know I'm performing safe injections every time. I think that ultrasound, particularly Clarius, is going to become the gold standard for filler treatments."

Dr. Steven F. Weiner - Facial Plastic Surgeon

World's First Ultra-High Frequency Scanner for Aestheticians



Clarius Cloud

Manage your exams and build reports anywhere, anytime with unlimited Clarius Cloud storage.



Easy To Use

Friendly AI-powered app for your Apple or Android automatically delivers the best imaging.



High-Definition

World's only 20 MHz wireless scanner for ultra-high resolution superficial imaging to 4 cm.



Specialized

AI-powered with specialized workflows for facial aesthetics.



Ultra-Portable

Ideal for small spaces, Clarius fits in your pocket. Take it anywhere!





Wireless

Steadfast connection. Total freedom. Without cords, cleaning and disinfection are simple and fast.




I trust Dr. Verner and I'm very happy that she's using ultrasound when she's treating my face. It makes me feel relaxed because I know I don't have to be concerned about something going wrong and I'm sure it's going to be perfect results."

Michal Paz - Patient of Dr. Ines Verner



Advanced Aesthetics Package

The Advanced Aesthetics Package offers more flexibility for users who need additional customizations for examinations and procedures in the realm of facial aesthetics, breast augmentation, BBL, hair follicle imaging, and more.



Dedicated Presets

Optimized presets for aesthetics, plastic surgery, and dermatology.



Vessel Depth AI

Automatically measure the shallowest flow found in real-time.



Split Screen

Add convenience by comparing left and right anatomy.

Additional Features

- DICOM
- Elastography
- Needle Enhance
- Pulsed Wave Doppler

Clarius Membership

The Aesthetics Advanced Software Package is Included with Your Clarius Membership

Membership Bundle

With your membership plan, you'll receive a discount on a scanner and the following:

- ✓ Support
- ✓ 60 Minutes of 1:1 Training
- ✓ Unlimited Clarius Cloud Storage
- ✓ Clarius Live
- ✓ Clarius Classroom
- ✓ Advanced Aesthetics Package

One-Time Purchase

You also have the option of paying for your scanner outright and buying option à la carte:

- ✓ Support
- ✓ 60 Minutes of 1:1 Training
- \$ Advanced Aesthetics Package

[REQUEST A QUOTE](#)

Clarius is on a mission to make accurate, easy-to-use, affordable ultrasound available to all medical professionals in every specialty. Our ultra-portable wireless scanners deliver best-in-class imaging with an easy-to-use app for iOS and Android that is powered by artificial intelligence and connected to the cloud. More than two million high-definition scans have been performed using Clarius wireless scanners, now available in over 90 countries worldwide.

Clarius Mobile Health

2985 Virtual Way #130,
Vancouver, BC, V5M 4X7
Canada
hello@clarius.com
Phone 1-778-800-9975

MKTG-00207-Rev 1