

ComfiDilator® 5-6 mm

D-52231-S




ComfiDilator® 5-6 mm

Comfi Dilator is designed to dilate the cervix in a gentle way. The dilator is 5 mm on one side and 6 mm on the other. The dilator is available in 2 sizes;

- Gentle
 - 2 sizes on the same instrument
 - Ergonomic thumb grip
 - Piece packaged
 - Sterile
 - Also available as 3-4 mm
- Melting temperature: 100 ° C
 - Not autoclavable
 - Latex free

Material: HD Polypropylene



ARTICLE NO.	PRODUCT	DEPT.PACKAGING	TRANSP.PACKAGING	TRANSP.PACKAGING
D52231-S	7 g	5 pcs	50 pcs	500 g

