



ComfiDilator® 3-4 mm

D-52230-S




ComfiDilator® 3-4 mm


Comfi Dilator is designed to dilate the cervix in a gentle way. The dilator is 3 mm on one side and 4 mm on the other. The dilator is available in 2 sizes;

- Gentle
 - 2 sizes on the same instrument
 - Ergonomic thumb grip
 - Piece packaged
 - Sterile
 - Also available as 5-6 mm
- Smälttemperatur: 90-100 °C
 - Ej autoklaverbar
 - Latexfri

Material: Polypropylene



ARTICLE NO.	PRODUCT	DEPT.PACKAGING	TRANSP.PACKAGING	TRANSP.PACKAGING
D52230-S	5 g	5 pcs	50 pcs	500 g



CETRO MEDICAL

CETRO MEDICAL AB, Nitgatan 11, 333 33 Smålandsstenar, Sweden, +46 371 330 30, www.cetromedical.se