



Specification

DNA/RNA-ExitusPlus™

A7089

Product Code:	A7089
Product Name:	DNA/RNA-ExitusPlus™
Specifications:	: • removes DNA and RNA contaminations from surfaces : • for cleaning of PCR work stations and equipment : • for cleaning of electrophoresis equipment, pipetts, reaction tubes etc. : • delivered as spray bottles or as refill bottles (RF) : ♦ The special features : • non-enzymatic degradation of DNA and RNA by catalytic and cooperative effects of the components : • All components of DNA-ExitusPlus™ are readily biologically degradable and not harmful or toxic for humans. : • Doesn't contain aggressive mineralic acids or alkaline substances : • No toxic fumes. Contains low concentration of alcohol. - Functionality (Strand breakage test): passes test - Efficacy (ExitusPlus™ activity test): passes test
WGK:	1
Storage:	RT protected from light
CS:	38221900

AppliChem GmbH

Ottoweg 4 • D-64291 Darmstadt • Phone +49 6151 9357 0 • Fax +49 6151 9357 11 • info.de@itwreagents.com • www.itwreagents.com
 CEO Joan Roget • Commerzbank Darmstadt • Bank 508 800 50 • Account 0186989900 IBAN DE24 5088 0050 0186 9899 00 • Swiftcode DRESDEFF508 • Finanzamt Darmstadt 007 225 00192 • Register court Darmstadt HRB Nr. 7340