



1



Insert smartphone
Tighten thumbscrew

2



Center over camera
Tighten thumbscrew

dermLite Universal Smartphone Adapter

3Gen



DermLite and 3Gen are registered trademarks of 3Gen, Inc. | 31521 Rancho Viejo Road | Suite 104 San Juan Capistrano, CA 92675 | U.S.A.

Tel. +1-949-481-6384 | Fax +1-949-240-7492

E-Mail: info@3geninc.com | www.derm-lite.com

©2017 by 3Gen, Inc. | All rights reserved.

REF A DLUA-1701

Patents/Trademarks: www.derm-lite.com/patents



M. Devices Group, Healthcare Education Centre,
The Church, Portland Street, Southport, PR8 1HU, UK.
Tel: +44 1704 544 944 | Fax: +44 1704 544 050



55 mm min.

98 mm max.