


Art. No.	
NAME OF PRODUCT	Alcohol swabs
DESCRIPTION OF PRODUCT	Alcohol swabs with 70% Isopropyl Alcohol Size: 60×30mm Remark: Specifications pls. refer to attachments I, II
INTENDED USE OF PRODUCT	For cleansing the skin

RAW MATERIAL		
CODE FOR THE ENCLOSED	RAW MATERIAL	ADDITIONAL INFORMATION
RAW MATERIAL SPECIFICATION	Non woven	Spunlace, white Weight: 50g/m ²
	Liquid	70% Isopropyl Alcohol
	Inner bag	Tinfoil paper
	Closure tape	PP transparent tape
	Inner box	Cardbord , white Weight: .350g/m ²
	Transport carton	Double Wall Corrugated Fibreboard ,white Weight: 600-700g/m ²

QUALITY SPECIFICATION OF THE PRODUCT			
CHARACTERISTIC	TARGET VALUE	QUALITY TOLERANCES	TEST METHOD
FINISHED SIZE	Length=60mm	+/- 2mm	Measuring
	Width=30mm	+/- 2mm	Measuring
FINISHED WEIGHT	Weight= 0.9g/pc	+/- 0.1g	Weighing

DESCRIPTION OF THE PACKAGING SYSTEM			
PACKING			
1pc/pouch, 100pcs/box, 50boxes/ctn			
INNER BAG SIZE (cm)			
L	W	T	
5	5	75mic	
BOX SIZE (cm)			
L	W	H	
10.3	5.5	5.2	
CARTON SIZE (cm)			
L	W	H	
29.5	22.5	28	
 The diagram illustrates the packaging process. It starts with a stack of individual pouches, labeled '×100pouches'. An arrow points to a box, which is labeled 'x50boxes'. Another arrow points down to a larger carton, also labeled 'x50boxes'.			

Laminate Alcohol swab

Attachment I

