

NEW



POLYCARE BIODEGRADABLE APRONS

CE certified, high quality single-use disposable aprons

As part of our commitment to the environment and sustainability, GV Health have introduced a new range of Polycare Biodegradable Aprons. Biodegradable aprons ensure speedier breakdown when discarded into landfill, where the majority of disposable polythene aprons currently are sent after use.

KEY FEATURES

Polycare biodegradable aprons are adjustable to fit all shapes and builds. Supplied on rolls for easy dispensing. Polycare biodegradable aprons are for use in all departments and medical specialisations including nursing, midwifery, infection control, domestic, paediatric and theatre.

- Standard quantity = 5 rolls of 200 aprons in biodegradable sleeves
- CE Certified
- Boxed in cardboard outers which can be recycled.
- 18 - 36 months to commence break down once sent to landfill

NEW POLYCARE RANGE

The aprons are just one of our new ranges we are bringing to market with sustainability at the forefront of its design.

DUTY OF CARE TO THE ENVIRONMENT

GV Health has taken important steps to reduce the amount of virgin polythene and significantly increase the amount of recycled material across our complete range of healthcare disposables.

• less polythene pollution of seas and seashores improves wildlife environments and ecosystems

• reduces the need for incineration, currently a major challenge in the UK due to limited number of sites that can recycle the large tonnage of plastic used

China was the largest importer of the UK's plastic waste, taking millions of tonnes of plastic a year. Since October 2017 China has ended all imports of UK plastic waste leaving the UK with a serious problem in how to process and recycle plastic. Currently the only option in many areas is landfill.

The UK Government is investing in new recycling projects but is a long way from being able to manage the quantity of plastic waste currently produced. 2042 is the current date for achieving the UK targets. GV Health, as a main approved supplier to the NHS is working with all partners to ensure the most sustainable way forward in terms of product and disposal NOW.



NHS SC Code	Direct Code	Description	Packing	Units per Roll/ Pack	Rolls/ Packs per case	Units per case	Cases per Pallet
BIODEGRADABLE DISPOSABLE APRONS							
BTB116	AWRB4616	Apron - Size 27 x 46 inches, 16mu White unisex Biodegradable disposable	Roll	200	5	1000	96
STANDARD APRONS							
BTB117	AYR4616	Apron - Size 27 x 46 inches, 16mu Yellow unisex disposable apron	Roll	200	5	1000	96
BTB118	AGR4616	Apron - Size 27 x 46 inches, 16mu Green unisex disposable apron	Roll	200	5	1000	96
BTB119	ABR4616	Apron - Size 27 x 46 inches, 16mu Blue unisex disposable apron	Roll	200	5	1000	96



Biodegradable



GV Health Ltd © T +44 (0)1920 463 098 F +44 (0)1920 484 664 E info@gvhealth.com
Outside UK: E export@gvhealth.com
www.gvhealth.com