

SEVERO

The best professional
pill grinder for hospitals,
nursing homes and home
care organizations.



SURGICAL HOUSE
Healthcare Consumables & Equipment

Key Features

- ✓ Fast and effective crushing of medication
- ✓ Avoiding exposure to medication and medication dust
- ✓ Small and portable design
- ✓ Easy to Clean
- ✓ Near silent operation
- ✓ Rechargeable battery
- ✓ Produced in the Netherlands

Products/Article codes

Severo Electric pill grinder (Art. code: SEVERO3.0), including

- Stainless steel grinding head
- Power supply

Recharger Disposables

- Severo cups (Art. code: MED001): 2100 cups per box
- Severo rondels (Art. code: MED003): 2100 cups per box
- Optional:

Severo cup lids (Art. code: MED002): 2000 lids per box

Grinding Specifications

Grinding type	Pill	Autonomous grinding cycles via combined compression and rotation.
grinding	Cup	Grinding of multiple standard pills and coated pills per cycle
volume	Grinding	30 ml
cycle time	Sound	8 seconds for regular pills
level	Battery last	52 dB (soft whisper level)
Charging	time	500 crushing's (= ca three days of continuous use)
Switches	Indicators	2 hours
Warranty		On/O
		Battery load and battery quality indicator
		2 years on the device, 1 year on the battery (not incl. user error)

Technical Specifications

Nominal	Voltage	12 Volts
Approximate	weight	ca 1,5 Kilograms
Dimensions	(LxBxH)	80 x 70 x 217 mm
Terminals	AC/DC	5,5 x 2.1 x 12 mm DC Connector for recharging
power supply	Battery	13 volts at 3 A
Battery	self-	Rechargeable sealed lithium-ion, 12 Volts 3 Ah
discharge	Fuse	ca 1% of capacity per month (25 c°)
		5A

Quality Assurance

Before shipment each device is thoroughly tested/inspected on:

- ✓ Scratches, dents, dirt or other exterior deformations
- ✓ Odd sounds during operation
- ✓ Functioning and grinding quality
- ✓ Battery capacity
- ✓ Availability of accessories

For further information:

