



ENCORE® Latex Micro

S U R G I C A L G L O V E S

Latex, powder-free surgical glove with exceptional dexterity and tactile sensitivity for surgical finesse

KEY FEATURES AND BENEFITS:

- 20% thinner on average than Ansell's standard latex powder-free surgical gloves with soft latex formulation for unsurpassed dexterity and tactile sensitivity
- Its brown color helps reduce glare of surgical lights and makes it suitable to be used as an underglove
- The DERMASHIELD™ Technology internal coating facilitates dry or damp donning and intra-operative changes

CAUTION: Safe use of these gloves by latex-sensitized individuals has not been established. This product contains natural rubber latex, which may cause allergic reactions.

Recommended For

- All surgeries, especially delicate procedures in ophthalmology, micro, and cardiovascular
- Damp donning
- Enhanced tactile sensitivity and dexterity, especially for delicate procedures
- Prevention of eye fatigue

FEATURED TECHNOLOGIES



ENCORE®
The Confidence to Perform

Dedicated to delivering optimized solutions engineered to provide a high degree of confidence to surgeons, O.R. nurses and surgical technicians seeking protection with a high level of dexterity, comfort and grip so they can perform at their best.



ENCORE[®] Latex Micro

SURGICAL GLOVES

PRODUCT INFORMATION

Material	Natural Rubber Latex Powder-Free Brown Micro
Powder Content	Beaded with SUREFIT [™] Technology Smooth with
Color	Micro-Textured Finish Polymer Coated with
Thickness	DERMASHIELD [™] Technology None Moderate
Cuff Style	
External Glove Surface	
Internal Glove Surface	
Allergy Prevention	
Grip Level	
Double Gloving Recommendation	Outer glove or underglove
Tested For Use With	
Chemotherapy Drugs Freedom from Holes (Inspection level I)	Yes, in accordance with ASTM D6978 (Not listed in the US FDA 510k)
	0.65 AQL
Product Standards	ASTM D3577, ASTM D7160, EN 16523-1, EN 374:2003, EN 420:2003 + A1:2009, EN 455 1-4, EN ISO 374-1:2016, EN ISO 374-5:2016, ISO 10282
Quality / Environmental Standards	EN 556, ISO 11137-Part 1, ISO 13485, ISO 14001, ISO 9001
Packaging	50 pairs per box; 4 boxes per carton/case; 200 pairs per carton/case
Shelf Life	3 years
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.
Vulcanization Chemical Accelerators	Zinc Dibutyldithiocarbamate (ZDBC), Zinc Diethyldithiocarbamate (ZDEC)

ORDERING INFORMATION

Size	5.5	6	6.5	7	7.5	8	8.5	9
Product Code	330104055	330104060	330104065	330104070	330104075	330104080	330104085	330104090

Contact your Ansell representative for ordering or more information.

North America
US Tel: 800 952 9916
CA Tel: 844 494 7854
Email: insidesalesus@ansell.com
Sandel US Tel: 866 764 3327
Sandel email: sandel-cs@ansell.com

Central & South America
Tel: +52(442) 296 20 50
Email: cslac@ansell.com
Brazil
Tel: +55-11-3356-3100
Email: luvas.medicas@ansell.com

Europe, Middle East & Africa
Tel: +32 (0) 2 528 74 00
Email: info@ansell.eu

Asia Pacific
Tel: +603 8310 6688
Email: apac.medical@ansell.com

India
Tel: +91 22 38124500
Email: info.india@ansell.com
Australia & New Zealand
Tel: +61 3 9270 7270
Email: protection@ap.ansell.com
China
Tel: +86(21) 3827 5005
Email: infochina@ansell.com

ansell.com

Ansell, and [®] and [™] are trademarks owned by Ansell Limited or one of its affiliates. US Patented and US and non-US Patents Pending:
www.ansell.com/patentmarking.
© 2021 Ansell Limited. All Rights Reserved.