

NEW ASPIRET

DESK-TYPE

A surgical suction unit indicated for minor surgical applications and post-operative therapy.

Thanks to its compact size it finds application in out-patient and in-patient care, elderly care and in private care.

The device is suitable for nasal, oral and endotracheal aspiration of bodily fluids (mucus or catarrh) from adults and children.

MAIN FEATURES



Non-stop operation

Approved for non-stop operation, no interruption to cool down is required.



Small size - high performance

Powerful engine, reaching quickly a max vacuum of 0,75 Bar.



High-quality materials

The chassis is made of VO(UL94) ABS durable plastic, flame retardant, corrosion resistant and rust-proof.



Long lasting accessories

Made of resistant, autoclavable materials to guarantee safety and hygiene over time.



Easy to handle

Tiny and lightweight, can be placed anywhere and carried with ease.



Best Seller configuration

NEW ASPIRET +1000 mL reusable jar
in Makrolon Polycarbonate (121°C)

Product code: RE 310001

OTHER CONFIGURATIONS

	1000 mL	2000 mL
REUSABLE SYSTEM (MAKROLON JARS)	RE 310001	RE 310001/01
DISPOSABLE SYSTEM (FLOVAC)	RE 310002	RE 310002/01

The warranty lapses immediately if the unit is used without all the required accessories or with non-original CA-MI accessories.



APPLICATIONS

HOME-CARE	GENERAELN STURGERY
RHOOME /- CWAARRID	GTERNAECRHAEL OSSU

TECHNICAL DATA

MOTOR	Oilless and maintenance-free piston pump
POWER FEEDING	220-230 V / 50 Hz or 120-127 V / 60 Hz
ISO 10079-1 CLASSIFICATION	High vacuum / Low flow
MAX VACUUM	-0,75 Bar -75 kPa -563 mmHg <i>adjustable value at sea level different altitudes may affect it</i>
MAX FREE AIR FLOW RATE <i>without accessories</i>	15 L/min
NOISE LEVEL	59,6 dB(A)
DUTY CYCLE	Non-stop operation
IP PROTECTION DEGREE	IP21
WEIGHT	2,20 Kg
DIMENSIONS	19 x 23 x 15 cm
WARRANTY	2 years
SHIPPING CARTON	4 units
MANUFACTURING COUNTRY	Italy



CONVENIENT HANDLE



INCLUDED ACCESSORIES

REUSABLE JAR OR DISPOSABLE SYSTEM

depending on configuration

REUSABLE SILICONE TUBES

Ø6 / 10 mm

BI-CONICAL CONNECTORS

ø 8 / 9 / 10 mm

ANTIBACTERIAL & HYDROPHOBIC FILTER

CH20 SUCTION CANULA

POWER CORD WITH SCHUKO PLUG

other plugs available upon request

OPTIONAL ACCESSORIES

Reusable Jars in Makrolon, APEC
Polycarbonate or Polysulfone

AVAILABLE SIZES:
1000 mL / 2000 mL

AVAILABLE SETS:

- ›Complete set
(Jar+cover with overflow system+ O-ring)
- ›Jar (without cover)
- ›Cover with overflow system
- › O-ring

For further details consult the Reusable
Collection System Table



OPERATING WITH:
REUSABLE COMPLETE TUBES' SET
3 silicone tubes (Ø 6/10mm, L: 140cm/ 13cm each)
+1 bi-conical connector Ø 8-9-10mm
+1 Antibacterial and hydrophobic filter




Disposable Flovac® System

AVAILABLE SIZES:
1000 mL / 2000 mL

AVAILABLE SETS:

- ›Reusable Flovac® Canister
- ›Disposable Flovac® Liner
- ›Male Connectors for Flovac® Liner
(5's pack)

For further details consult the Disposable
Collection System Table



OPERATING WITH:
FLOVAC COMPLETE TUBES' SET
2 tubes (Ø 6 / 10mm, L: 140cm / L: 26cm)
+1 bi-conical connector (Ø 8 / 9 / 10mm)



Rolls of silicone tube

AVAILABLE SIZES:

- ›Ø 6 / 10mm, Length = 1m
- ›Ø 6 / 10mm, Length = 10m
- ›Ø 6 / 10mm, Length = 50m



Cannulas

AVAILABLE MODELS

- ›Aspiration Probe CH20 (set of 10's)
- ›Yankauer Handle Flat Tip With Hole
- ›Yankauer Handle Crown Tip With Hole
- ›Yankauer Tube L= 180cm



Stands

Stand on 5 Castor Ø 50mm with brakes



Carrying Bags

AVAILABLE SIZES:
1000 mL / 2000 mL

