

Technical Data Sheet

Article-No.: **01035**
Description: **MED COMFORT**
Latex finger-cots
white, non sterile, powder free

MED-COMFORT®
Latex-Fingerlinge puderfrei
Latex finger cots powder free

100 Stück/Pieces
Größe/Size
2 ☐ Art.-No. 01035-2
3 ☐ Art.-No. 01035-3
4 ☐ Art.-No. 01035-4
5 ☐ Art.-No. 01035-5

gerollt nicht steril
roll-up non-sterile
© 2019 Günter + Berndtson AG
21423 Winsen (Luhe) - Germany
Wichtiger Hinweis: Dieses Produkt enthält Naturkautschuklatex, das allergische Reaktionen auslösen kann. Personen, die allergisch gegen Naturkautschuk sind, sollten dieses Produkt nicht verwenden. Lesen Sie bitte die Gebrauchsanweisung. Für medizinische Risiken: Für Latex-Nutzung.

PRAXISDIENST
Medical Supplies since 1953
Order here!

THICKNESS
cuff
0,11 mm

PRODUCT DESCRIPTION							
material	☒ Latex	☐ Nitrile	☐ Vinyl	☐ Vinyl-Nitrile-mixture	☐ Polyethylen (PE)	☐ TPE	☐ cotton
colour	☒ white	☐ blue	☐ black	☐ mint	☐ purple	☐ mix	☐ bordeaux
characteristics	☐ pre-powdered	☒ powder free	☐ sterile	☒ non sterile	☒ ambidextrous	☐ fits hand-specific	
surface	☐ fully textured	☐ finger textured	☒ non textured	☐ embossed	☐ outside chlorinated, coated with Aloe Vera inside		

SIZES							
	XS (1)	S (2)	M (3)	L (4)	XL (5)	XXL (6)	XXXL (7)
width	-	24 ± 2 mm	27 ± 2 mm	30 ± 2 mm	34,5 ± 2 mm	-	-
length	-	70 ± 3 mm	70 ± 3 mm	70 ± 3 mm	70 ± 3 mm	-	-

REGULATORY AFFAIRS							
PPE-Regulation (EU) 2016/425	☐ Category I	☐ Category II	☐ Category III	☒ no PPE-article			
MD-Regulation (EU) 2017/745	☐ Class I	☐ Class II	☐ Class III	☐ sterile	☐ measuring function	☒ no medical device	
Food Contact (EG) 1935/2004	☐ acidic foods	☐ aqueous foods	☐ fatty foods	☐ alcoholic foods	☐ dry foods	☒ not approved for food-contact	

STANDARDISATION							
EN 388:2016 Mechanical Risks	abrasion resistance	blade cut resistance Coupe-Test	tear resistance	puncture resistance	blade cut resistance TDM-Test	impact test	
Level	not applicable						

EN 374-1:2016 Chemical Risks	chemical	code letter	level	permeation	degradation	
	not applicable					
EN 374-4:2013 Degradation						
EN 374-5:2016 microorganism tightness	not applicable					

EN 420:2010 protective gloves	not applicable					
-------------------------------	----------------	--	--	--	--	--

EN 455 medical gloves for single use	not applicable					
EN 455-1:2000 freedom from holes	not applicable					

Technical Data Sheet

Article-No.: **01035**
Description: **MED COMFORT**
Latex finger-cots
white, non sterile, powder free



LOGISTIC DATA SUBPACKING		
generell information		
material	carton	
pieces per subpacking	100	
EAN subpacking size XS	-	
EAN subpacking size S	-	
EAN subpacking size M	-	
EAN subpacking size L	-	
EAN subpacking size XL	-	
EAN subpacking size XXL	-	
EAN subpacking size XXXL	-	
PZN subpacking size XS	-	
PZN subpacking size S	-	
PZN subpacking size M	-	
PZN subpacking size L	-	
PZN subpacking size XL	-	
PZN subpacking size XXL	-	
PZN subpacking size XXXL	-	
measures & size		
length	240 mm	
width	120 mm	
height	70 mm	
weights		
size	net weight	gross weight
XS	-	-
S	42 g	52 g
M	47 g	57 g
L	52 g	62 g
XL	57 g	67 g
XXL	-	-
XXXL	-	-

LOGISTIC DATA PALETTE		
general information		
kind of palett		euro-palette
measures & size		
cartons per layer		4
layers per palette		5
height of the palette		158 cm
weights		
size	net weight	gross weight
XS	-	-
S	322 kg	347 kg
M	352 kg	377 kg
L	382 kg	407 kg
XL	412 kg	437 kg
XXL	-	-
XXXL	-	-

LOGISTIC DATA OUTER PACKING		
generell information		
material	carton	
subpackings per outer packing	300	
EAN outer packing size XS	-	
EAN outer packing size S	4044941000770	
EAN outer packing size M	4044941000787	
EAN outer packing size L	4044941000794	
EAN outer packing size XL	4044941000800	
EAN outer packing size XXL	-	
EAN outer packing size XXXL	-	
PZN outer packing size XS	-	
PZN outer packing size S	-	
PZN outer packing size M	-	
PZN outer packing size L	-	
PZN outer packing size XL	-	
PZN outer packing size XXL	-	
PZN outer packing size XXXL	-	
measures & size		
length	490 mm	
width	400 mm	
height	285 mm	
weights		
size	net weight	gross weight
XS	-	-
S	15.600 g	16.100 g
M	17.100 g	17.600 g
L	18.600 g	19.100 g
XL	20.100 g	20.600 g
XXL	-	-
XXXL	-	-

Rev.-Nr.: 5
Datum 02. 03.2022

changes and errors excepted