

1

Otoskope
Otosscopes



LUXAMED®

Medizintechnik Made in Germany



2.5 V

A1.416.114	schwarz black	■
A1.416.214	blau blue	■
A1.416.314	grau grey	■
A1.416.914	weiß white	□



LuxaScope Auris LED 2.5 V

- innovatives Otoskop im Taschenformat
- Beleuchtung: LED-Ring mit 6 LEDs (konzentrisch), keine Fiber-Optik, kein eingeschränktes Sichtfeld, Geschützt durch DE202013104278U1
- wartungsfrei, da kein Lampenwechsel notwendig
- dimmbar
- Abschaltautomatik nach 3 Minuten
- Lebensdauer LEDs: ca. 100.000 Stunden durch Konstantstromansteuerung
- Beleuchtungsstärke: ca. 10.000 Lux
- Farbtemperatur: ca. 4.200 K
- optimierte Bauform des Kopfes
- schwenkbare 3-fach-Vergrößerungslupe mit Rastfunktion
- Gebläseanschluss für pneumatischen Test
- 2-Komponenten-Griff (Aluminium / glasfaserverstärkter Kunststoff), kombinierbar auch mit LuxaScope LED Dermatoskop- und Ophthalmoskop Kopf
- elektronischer ON/OFF Druckschalter
- Stromversorgung: 2 x AAA Alkaline Batterien, inklusive
- inkl. 20 Einweg-Ohrtrichter, grau, je 10 x 2.5 mm und 10 x 4.0 mm

LuxaScope Auris LED 2.5 V

- innovative otoscope in pocket size
- illumination: LED-Ring with 6 LEDs (concentric), no fibre optic, field of vision is not restricted, Protected by utility patent DE202013104278U1
- maintenance-free, because no lamp replacement necessary
- can be dimmed
- automatic switchoff after 3 minutes
- life expectancy of LEDs: approx. 100.000 hours by constant current drive
- intensity of illumination: approx. 10.000 Lux
- colour temperature: approx. 4.200 K
- optimized design of head
- pivoting 3-fold magnifying glass with locking function
- insufflation port for pneumatic test
- two-component handle (aluminium / fibre-glass reinforced), also combinable with LuxaScope LED Dermatoscope- and Ophthalmoscope head
- electronic ON/OFF push-button switch
- power supply: 2x AAA Alkaline batteries, inclusive
- incl. 20 disposable ear funnels, grey, each 10 x 2.5 mm and 10 x 4.0 mm



A1.416.934	weiß / weißblau	white / ice-blue	■
A1.416.944	weiß / rapsgelb	white / colza yellow	■
A1.416.954	weiß / jadegrün	white / jade-green	■
A1.416.964	weiß / zartflieger	white / delicate lilac	■
A1.416.974	weiß / magenta	white / magenta	■