

EasyOne FlowTube

Breathing mouthpiece for Infection Control



n d d **EasyOne FlowTube** mouthpieces avoid cross-contamination between patients and reduce cleaning of the device to an absolute minimum.

The **EasyOne FlowTube** mouthpieces have been engineered to fit perfectly into n d d's device to guarantee reliable and high quality measurement results.

All parts exposed to patient's breath are **single-use and disposable**

No expiration date

Technical Data

Material

EasyOne FlowTube: Polypropylene (PP)
Bag: Polyethylene (PE)

Transport and storage information

Temperature: -25 to 70°C
Rel. Humidity: 5 to 90%, non-condensing
Atmospheric Pressure: 400 to 1060 hPa

Dimension (single EasyOne FlowTube)

21 x 38 x 150 mm (H x W x D), 8 g

Compatible with

EasyOne Air

Order Information

Part number:

5050-50
5050-200
5050-500

Product:

EasyOne FlowTube, individually wrapped, dispenser box of 50 pcs
EasyOne FlowTube, individually wrapped, box of 200 pcs (dispenser box included)
EasyOne FlowTube, individually wrapped, box of 500 pcs (without dispenser box)