

# EasyOne FlowTube

Breathing mouthpiece for Infection Control

**n d d**  
Medical Technologies



n d d **EasyOne FlowTube** mouthpieces avoid cross-contamination between patients and reduce cleaning of the device to an absolute minimum.

The **EasyOne FlowTube** mouthpieces have been engineered to fit perfectly into n d d's device to guarantee reliable and high quality measurement results.

All parts exposed to patient's breath are **single-use and disposable**

No expiration date

## Technical Data

### Material

EasyOne FlowTube: Polypropylene (PP)  
Bag: Polyethylene (PE)

### Transport and storage information

Temperature: -25 to 70°C  
Rel. Humidity: 5 to 90%, non-condensing  
Atmospheric Pressure: 400 to 1060 hPa

### Dimension (single EasyOne FlowTube)

21 x 38 x 150 mm (H x W x D), 8 g

### Compatible with

EasyOne Air

## Order Information

### Part number:

5050-50  
5050-200  
5050-500

### Product:

EasyOne FlowTube, individually wrapped, dispenser box of 50 pcs  
EasyOne FlowTube, individually wrapped, box of 200 pcs (dispenser box included)  
EasyOne FlowTube, individually wrapped, box of 500 pcs (without dispenser box)