

Spirette™

Breathing mouthpiece for Infection Control



Medical Technologies



Order here!



The ndd Spirette mouthpieces have been engineered to **fit perfectly into ndd devices to guarantee reliable and high quality measurement results.**

ndd Spirette mouthpieces **prioritise hygiene and safe re-use of ndd devices** by avoiding cross-contamination between patients.

All parts exposed to patient's breath are **single-use and disposable**

No expiration date

Technical Data

Material

Spirette: linear low density polyethylene (LLDPE)
Bag: Polyethylene (PE)

Transport and storage information

Temperature: -25 to 70°C
Rel. Humidity: 5 to 90%, non-condensing
Atmospheric Pressure: 400 to 1060 hPa

Dimension (single Spirette)

25 x 38 x 140 mm (H x W x D), 8 g

Compatible with

EasyOne, Easy on-PC, EasyOne Pro and EasyOne Pro LAB

Order Information

Part number:

2050-1

2050-5

2050-10

Product:

Spirette, individually wrapped, dispenser box of 50 pcs.

Spirette, individually wrapped, box of 200 pcs. (dispenser box included) wrapped, bag of 10 pcs.

Spirette, individually wrapped, box of 500 pcs. (without dispenser box)