

Pen options

NovoFine® needles fit all Novo Nordisk, and most diabetes pens.



FlexTouch®



NovoPen® 5



FlexPen®



NovoPen® 4



NovoPen Echo®

NovoFine® needles also fit most diabetes pens, such as:

- | | | |
|------------------|--------------------|-----------------------|
| KwikPen® | • Autopen® 24 | • BerliPen® 302 |
| SoloStar® | • Omnicar® Pen 31 | • HumaPen® Luxura™ |
| Liprolog® Pen | • OptiPen® Pro 1 | • HumaPen® Luxura™ HD |
| Huminsulin® Pen | • Ypsopen® | • HumaPen® MEMOIR™ |
| Exenatide Pen | • BerliPen® Areo | • OptiClik® |
| Autopen® Classic | • BerliPen® Areo 2 | • TactiPen® |

Brands are trademarks of their respective owners.

The NovoFine® range

Available in packs of 100

Length	Diameter
4 mm	32G*
6 mm	32G*
6 mm	31G*
8 mm	30G*



8 mm

30G* Autocover®

* The G stands for gauge and is a measure of needle thickness; the higher the number, the thinner the needle.