

SURGICAL MESH

Product

Surgical mesh is a medical device that is used to provide additional support to weakened or damaged tissue.

Indications

Reinforcement of the abdominal wall hernia / eventration via celioscopy or laparotomy

Attributes

- ✓ No shrinkage
- ✓ No coating
- ✓ Ethylene Oxide sterilisation
- ✓ Non-absorbable

Benefits

- ✓ Extreme dimensional stability
- ✓ Very thin mesh structure
- ✓ Ideal porosity – well tolerated
- ✓ Excellent transparency
- ✓ Good handling - easy to cut

Technical composition

Monofilament polypropylene knitted into an elastic, durable mesh

Square & Rectangular Mesh	Size	Order no
standard mesh	5 cm x 10 cm	ZMS0510
5 pieces / pack	6 cm x 11 cm	ZMS0611
	7,5 cm x 15 cm	ZMS07515
	10 cm x 15 cm	ZMS1015
	15 cm x 15 cm	ZMS1515
	30 cm x 30 cm	ZMS3030



*For quality
&
safety*

smi
www.sutures.be