

	NATURAL RUBBER LATEX (NR)						
BRAND NAME	sempermed® supreme		sempermed® supreme plus		sempermed® supreme green		
 quick don™  quick don™  quick don™							
PRODUCT BENEFITS	A PREMIUM ALL-ROUNDER YOU CAN TRUST ■ Extra-supple: great comfort & tactile sensitivity ■ Fast & easy to don: Quick Don™ inner coating ■ Anatomical shape: reduces strain for wearer ■ Skin-friendly: washed to minimise skin irritants		OPTIMAL GRIP IN EVERY SITUATION ■ Micro-rough surface: for dry & damp applications ■ Fast & easy to don: Quick Don™ inner coating ■ Anatomical shape: reduces strain for wearer ■ Skin-friendly: washed to minimise skin irritants		DOUBLE PROTECTION THAT YOU CAN TRUST ■ Clear detection of perforations: when double-donning with a lighter coloured outer glove ■ Tactile sensitivity thanks to optimised wall thickness ■ Ideal inner glove: designed for double-donning ■ As underglove for sempermed® supreme or supreme plus		
TYPE	Sterile surgical glove for single use, powder free, with synthetic lining						
MATERIAL	Natural rubber latex (NR)		Natural rubber latex (NR)		Natural rubber latex (NR)		
COLOUR	Natural white		Natural white		Green		
GLOVE SHAPE	Fully anatomical with curved fingers and rolled rim						
SIZE/OVERALL LENGTH as per EN 455-2			5 ½, 6 and 6 ½ median 270 mm 7, 7 ½ and 8 median 280 mm 8 ½ and 9 median 285 mm				
WALL THICKNESS measured in single layer in the palm area	0.19–0.23 mm		0.19–0.23 mm		0.17–0.20 mm		
BARRIER PERFORMANCE	AQL 0.65						
FORCE AT BREAK as per EN 455-2	≥ 12 N, exceeding EN 455-2 (≥ 12 N instead of 9 N)						
DURABILITY in original packaging if stored as per DIN 7716, ISO 2230	3 years						
STERILISATION	Radiation STERILE R						
PACKAGING	Left and right glove with turned up cuff in inner pouch, ozone-tight, sealed in flat medical peel pack. In dispenser carton with sterilisation indicator: 50 pairs. In transport carton with sterilisation indicator: 300 pairs.						
PURPOSE according to MP Regulation 93/42/EWG / according to PPE Regulation (EU) 2016/425	Medical device class IIa Single-use protective glove – PPE category III (protection against chemical substances for limited time) ²⁾						
ARTICLE NUMBERS	EU		INT		EU		
	Size 5 ½	822751521	822751525	Size 5 ½	822851521	Size 5 ½	823756521
	Size 6	822751601	822751605	Size 6	822851601	Size 6	823756601
	Size 6 ½	822751621	822751625	Size 6 ½	822851621	Size 6 ½	823756621
	Size 7	822751701	822751705	Size 7	822851701	Size 7	823756701
	Size 7 ½	822751721	822751725	Size 7 ½	822851721	Size 7 ½	823756721
	Size 8	822751801	822751805	Size 8	822851801	Size 8	823756801
	Size 8 ½	822751821	822751825	Size 8 ½	822851821	Size 8 ½	823756821
	Size 9	822751901	822751905	Size 9	822851901	Size 9	823756901

¹⁾ Crepy et.al. (2018), Accelerator-free gloves as alternatives in cases of glove allergy in healthcare workers, in Contact Dermatitis. 2018 Jan;78(1):28-32. doi: 10.1111/cod.12860. Epub 2017 Jul 27. <https://www.ncbi.nlm.nih.gov/pubmed/28748553>

²⁾ Please note the Sempermed Chemical Resistance list, available on <https://www.sempermed.com/knowledge-area/chemical-resistance/chemical-resistance-sempermed-brands/>