

# MC-980MA PLUS

## Multi-Frequency Segmental Body Composition Analyser with Sarcopenia Indicator

- Full segmental body composition analysis in 30 seconds.
- Integrated Microsoft Windows Operating System allows simple link ability to other compatible devices.
- 6 frequencies (1kHz, 5kHz, 50kHz, 250kHz, 500kHz, 1000kHz) allowing intra and extra cellular water.
- NEW Skeletal Muscle Mass Index for in-depth sarcopenis assessment.
- Data can be printed onto consultation sheet using any printers, can also use together with optional accessories including barcode scanner and NFC card reader.
- Once registered, all readings will be saved into in-built database for record and easy review.
- In-built software in Asian and European languages.
- 300kg weighing capacity with wide, low-profile platform.
- Easy to use touch screen display allows free standing use.
- Normal BMI and Fat Range can be set by client, desirable range will be adjusted automatically.
- Calibrated up to 300,000 uses with automatic calibration after each measurement.

