

PLASTIC MEASURING JUG | 500 ML



DESCRIPTION

Durable measuring jug made of transparent polypropylene, used for precise liquid measurement. The easy-to-read scale is durably imprinted. You can order the measuring jug in three different sizes.

PRODUCT DETAILS

- Plastic measuring jug
- Material: polypropylene (PP)
- Imprinted scale
- With spout and handle
- Stackable
- Available in the following sizes: 500 ml, 1,000 ml, 2,000 ml
- Delivered quantity: 1 piece

Overview

Article Number:
T135857

Category:
Lab Equipment

EAN:
4260547369367

Teqler is a trademark of NetMed S.à.r.l.

 NetMed S.à.r.l.
8, Am Scheerleck, 6868 Wecker, Luxembourg

 +352 267149 09  +352 267149 19  info@netmed.lu

Social Media

