

ComfiSound Sterile

D-52208-S




ComfiSound Sterile

ComfiSound is a flexible uterus probe with centimeter scale for accurate measurement of uterus. The probe also has a removable marker for easier measurement.

- Flexible
 - Latex free
 - Easy grip
 - Individually packaged
- Melting temperature: 100 ° C
 - Not autoclavable
 - Latex free

Material: HD polyethylene, LD polyethylene



ARTICLE NO.	PRODUCT	DEPT.PACKAGING	TRANSP.PACKAGING	TRANSP.PACKAGING
D52208-S	10,5 g	10 pcs	50 pcs	100 g

